



# Naída CI Operation Guide


This guide provides instructions on the essential operations to use your Naída CI sound processor.



For more Naida CI information download the MyNaida CI app on iTunes or Google Play or visit [AdvancedBionics.com.au](http://AdvancedBionics.com.au).

## Turn On the Naida CI

Slide a fully charged battery onto the Naída CI sound processor.



1 Hold the Naída CI sound processor in one hand. Rotate the cable upward to avoid bumping it upon battery placement.



2 Align the battery with the Naída CI housing side of the connector.



3 Slide the battery onto the Naída CI until it clicks into place.

When the battery is engaged the LED light will flash ORANGE to indicate battery status.\* Three to four ORANGE blinks indicate the battery is sufficiently charged to power the Naída CI sound processor. The ORANGE blinks will be followed by one quick GREEN blink to indicate the Naída CI is in program one.

\*Battery status LEDs not available with disposable battery cartridge.

## Turn Off the Naida CI

Remove the battery by sliding it towards the cable side of the Naida CI.

## Volume Control

If volume control is enabled, the loudness level can be adjusted. Pressing the top of the button increases the volume, while pressing the bottom of the button decreases the volume.

## Program Control

If volume control is enabled, the loudness level can be adjusted. Pressing the top of the button increases the volume, while pressing the bottom of the button decreases the volume.



## Batteries

There are four rechargeable battery options for the Naida CI.



## How to Charge Rechargeable Batteries

Slide the rechargeable battery onto the slide tracks on the charger until it clicks into place. The charging slot will light up green when charging and turn blue when fully charged.

