ADVANCED
BIONICS

POWERFUL CONNECTIONS

Naída™ CI M90 Sound Processor

CI-5293

instructions for use

English, Français, Español, Deutsch, Nederlands, Italiano

Contents

SYMBOLS TABLE	4
PURPOSE AND SCOPE OF IFU	5
LIMITATIONS AND CONTRAINDICATIONS	5
PRECAUTIONS, CAUTIONS, AND WARNINGS	6
UNDESIRABLE SIDE EFFECTS	11
PRODUCT DESCRIPTION AND PERFORMANCE CHARACTERISTICS	11
INDICATIONS FOR USE	12
INTENDED PURPOSE	13
<i>Intended Users</i>	14
<i>Intended Use Environment</i>	14
COMPATIBILITY	15
OPERATING INSTRUCTIONS	16
<i>Powering Naída CI M90</i>	16
<i>Attaching the Battery</i>	17
<i>Removing the Battery</i>	18
<i>Approved Naída CI M90 Power Supplies</i>	19
<i>Connecting the earhook to the sound processor</i>	20
<i>Removing the earhook</i>	22
<i>Attaching the headpiece cable to the sound processor</i>	24
<i>Removing the Slim HP cable from the sound processor</i>	25
<i>Naída CI M90 LEDs</i>	25
<i>Naída CI M90 Internal Alarms</i>	28
<i>Using the Naída CI M90 Multi function Button</i>	29

Connectivity Overview	30
<i>Managing Bluetooth Phone Calls with the Naída CI M90</i>	31
<i>Streaming Audio via Bluetooth with the Naída CI M90</i>	32
<i>Pairing the Naída CI M90 to Phonak Accessories</i>	32
<i>Switching between Multiple Audio Sources with the Naída CI M90</i>	33
<i>Placing the Naída CI M90 in Flight Mode</i>	33
CARE AND MAINTENANCE	34
<i>Storing the Naída CI M90</i>	34
<i>Recommended Operating and Storage Temperature, Humidity and Pressure Ranges</i>	34
<i>Cleaning and Maintenance</i>	34
<i>Expected Life</i>	35
<i>Ingress Protection (IP) Ratings</i>	35
DISPOSAL AND HAZARDOUS MATERIALS	36
AVAILABLE PRODUCT OPTIONS AND MODEL NUMBERS	36
CLINICAL BENEFITS	37
COUNSELLING AND FITTING RECOMMENDATIONS	43
SUMMARY OF SAFETY AND CLINICAL PERFORMANCE	46
GUIDELINES FOR IT SECURITY	46
<i>Intended Use Environment</i>	46
<i>System and Network Diagram</i>	47
<i>System Interfaces</i>	47
<i>Cybersecurity Configuration</i>	49
<i>Integrated Cybersecurity Functions</i>	49

<i>System, Configuration, and Data Backup/Restore</i>	50
<i>Cybersecurity Incident Response</i>	50
<i>Patches and Updates</i>	50
<i>Available Training</i>	50
<i>End of Cybersecurity Support</i>	51
<i>GUIDANCE AND MANUFACTURER'S DECLARATION</i>	51
<i>Electromagnetic Emissions</i>	51
<i>Electromagnetic Immunity</i>	52
<i>Separation Distances between RF Communications Equipment and Naída CI M90</i>	54
<i>Quality of Service</i>	55
<i>Wireless Safety</i>	56

SYMBOLS TABLE

Labeling symbols and their meanings

	European Community Mark of Conformity. Authorized to affix the CE Mark in 2021
	Date of Manufacture
	Manufacturer
	Model number
	Serial number
	European authorized representative
	Unique Device Identifier
	Medical Device
	Type of Protection: BF
	Suitable temperature range for transport and storage
	See Instructions for Use

	Fragile
	Suitable relative humidity range for transport and storage
	Caution (refer to “Cautions & Warnings” section for further information)
	Dispose of in accordance with applicable national and local regulations
	MR Unsafe
IPxx	Ingress Protection Rating

PURPOSE AND SCOPE OF IFU

This IFU is designed to help cochlear implant recipients, and their caregivers if applicable, to understand the use and care of the Naída CI M90 sound processor. It is also meant to be used by hearing care professionals trained in the fitting of cochlear implants to counsel the recipients on the use of the sound processors and to troubleshoot sound processor problems.

Limitations and Contraindications

As the Naída CI M90 sound processor is a component of the HiResolution Bionic Ear System, the following contraindications stated for the HiResolution Bionic Ear System are applicable: Deafness due to lesions of the acoustic

nerve or central auditory pathway; cochlear ossification that prevents electrode insertion; absence of cochlear development.

Precautions, Cautions, and Warnings

- *This device should be used only by the individual for whom it is prescribed.*
- **CHOKING HAZARD:** *Contains small parts that pose a hazard of inhalation, choking, or ingestion. Using your sound processor and accessories contradictory to their intended use (e.g., mouthing, chewing) may cause bodily harm. If any parts are swallowed or inhaled, consult a physician or hospital immediately.*
- *Ensure appropriate supervision when child is wearing the Naída CI M90 sound processor and accessories. Do not allow children to play with or leave them unattended with the sound processor or any of its accessories.*
- *Power supplies and battery chargers should be operated in an open area to ensure adequate airflow. While no injury cases have resulted, components may become hot during normal use or a fault condition. If the device's temperature results in discomfort or pain when touched, disconnect the power source and contact your local Advanced Bionics representative.*
- *Do not use any other power supply with the sound processor or accessories unless it is supplied by Advanced Bionics. If needed, please call Advanced Bionics for a power supply replacement.*

- *Do not use accessories when they are plugged into power sources such as wall outlets, or other power sources that are USB compatible (e.g., laptops).*
- *Do not allow leaking battery fluid to come into contact with skin, mouth, or eyes.*
- *Do not expose batteries to heat (e.g., do not store in direct sunlight or in a hot car).*
- *Remove batteries from your sound processor when they are drained to prevent damage from possible leaking.*
- *Do not dispose of batteries in fire.*
- *High levels of static electricity may damage electronic components of your sound processor or implant. Care should be taken to avoid exposure of the system to static electricity.*
- *Do not attempt to pair your sound processor to any accessories while driving or operating heavy machinery.*
- *Pair the sound processor to Bluetooth devices only in a secure environment.*
- *Unwanted paired Bluetooth devices can be overwritten by pairing with authorized devices or can be deleted during fitting by your cochlear implant professional.*
- *Do not stream audio to your sound processor while driving or operating heavy machinery.*
- *Use your sound processor for hands-free phone use only where permitted by law and only when you will not be distracted from the safe operation of your motor vehicle.*

- *Attenuating ambient sounds or adjusting environmental balance while streaming audio (either with the multifunction button or with the mobile application) may compromise situational awareness.*
- *Do not use your Advanced Bionics accessories when instructed not to use wireless electronic devices, such as on airplanes.*
- *Remove your sound processor and headpiece before entering a room where an MRI scanner is located, and prior to having any surgeries that involve electrocautery and diathermy.*
- *Remove your sound processor and consult your cochlear implant professional if uncomfortable sounds are heard or in case of discomfort, pain or skin irritation.*
- *If the sound processor or accessories become unusually hot, or warm, discontinue use immediately and contact Advanced Bionics or your cochlear implant professional.*
- *The large Slim HP Color Cap is for use by recipients older than 3 years of age, as this color cap and the Slim HP magnet may become displaced upon dropping.*
- *It is important to have the correct magnet strength so you do not experience discomfort or retention issues. If magnet strength is insufficient, the headpiece may fall off more than is acceptable. If magnet strength is excessive, you may experience irritation or discomfort. Consult your cochlear implant professional if there are any concerns regarding magnet strength. If deemed appropriate, an audiologist may change the strength of the magnet in the headpiece. Do not change magnet strength unless under the direction*

of a cochlear implant professional. If you experience any redness, irritation, or discomfort, discontinue use of the headpiece immediately and contact a cochlear implant professional. See the headpiece instructions for use for additional information regarding adjustment of headpiece magnet strength.

- Store additional headpieces away from items with magnetic strips (e.g., credit cards, hotel room keys, etc.) as they may de-magnetize cards.*
- The digitally-coded, inductive transmission technology used in this device is extremely reliable and experiences virtually no interference from other devices. It should be noted, however, that when operating the device near a computer terminal or other strong electromagnetic field (e.g. RFID system, metal detectors, electromagnetic anti-theft systems), it may be necessary to be at least 24" (60 cm) away to ensure proper operation. If the Naída CI M90 sound processor does not respond to the implant device because of an unusual field disturbance, move away from the disturbing field.*
- The Naída CI M90 sound processors and accessories should be used in accordance with the electromagnetic compatibility information provided in the Guidance and Manufacturer's Declaration section of this document.*
- Although electromagnetic emissions of your Naída CI M90 sound processor have been confirmed within safe limits, some other devices could be sensitive to such emissions. If you notice other devices do not behave as expected when your sound processor is nearby, separate the devices.*

- *Portable and mobile RF communications equipment, including radios and cellular phones, may affect sound quality of the Naída CI M90 sound processor and accessories; however, there is no safety hazard associated with such equipment.*
- *Do not expose any part of the Naída CI M90 sound processors or accessories to extreme heat, such as an oven, microwave, or hair dryer.*
- *AutoSense OS™ 3.0 operating system, Speech in Loud Noise, Fixed Directional mode, UltraZoom, UltraZoom + SNR Boost, and Speech in 360 programs may dampen sounds that are not in front of the recipient.*
- *Use of WindBlock, EchoBlock, SoundRelax, NoiseBlock, and/or WhistleBlock may affect the quality of sound.*
- *Do not use any programs other than an off-ear program when utilizing any off-ear wearing configuration (e.g., in a clip, in the M Waterproof Battery).*
- ***NOTE:** In the United States, Fixed Directional mode, StereoZoom, UltraZoom, SoundRelax, NoiseBlock, WindBlock, and EchoBlock are approved for use in pediatric recipients 6 years and above who are 1) able to complete objective speech perception testing in order to determine speech performance and 2) able to report a preference for different coding strategies or features.*
- *In the event that you experience any issues with your product, please contact your cochlear implant professional or the manufacturer. Do not attempt to*

service or modify the Naída CI or its accessories. Doing so may compromise system performance and will void the manufacturer's warranty. Products should be serviced only at Advanced Bionics and damaged products should be returned to Advanced Bionics.

Undesirable Side Effects

Undesirable side effects of your Naída CI M90 system may include skin irritation and discomfort from pressure on the ear, device overheating, or overly loud sounds. If any undesirable side effect is encountered, please remove your sound processor and consult your cochlear implant professional.

Product Description and Performance Characteristics

The Naída CI M90 Sound Processor is a behind-the-ear (BTE) sound processor for use with an Advanced Bionics cochlear implant. The Naída CI M90 is comprised of the components pictured below.

The Naída CI M90 sound processor features direct connectivity to wireless accessories and Bluetooth®* technology enabled devices. The Naída CI M90 offers multiple power options and wearing options for the user.

The system is designed to provide useful hearing to individuals with severe-to-profound hearing loss.

INDICATIONS FOR USE

The Naída CI M90 sound processor is an external component of the HiResolution Bionic Ear System, which is intended to restore a level of auditory sensation to individuals with severe-to-profound sensorineural hearing loss via electrical stimulation of the auditory nerve.

Adults

- *18 years of age or older.*
- *Severe-to-profound bilateral sensorineural hearing loss or severe-to-profound unilateral hearing loss.*
- *Post-lingual onset of severe or profound hearing loss.*
- *Limited benefit from appropriately fitted hearing aids, defined as scoring 50% or less on a test of open-set sentence recognition (HINT Sentences).*

Children

- *12 months through 17 years of age.*
- *Severe-to-profound bilateral sensorineural deafness or severe-to-profound unilateral hearing loss.*

*The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.

- *Use of appropriately fitted hearing aids for at least 6 months in children 2 through 17 years of age, or at least 3 months in children 12 through 23 months of age. The minimum duration of hearing aid use is waived if x-rays indicate ossification of the cochlea.*
- *Little or no benefit from appropriately fitted hearing aids. In younger children (< 4 years of age), lack of benefit is defined as a failure to reach developmentally appropriate auditory milestones (such as spontaneous response to name in quiet or to environmental sounds) measured using the Infant-Toddler Meaningful Auditory Integration Scale or Meaningful Auditory Integration Scale or $\leq 20\%$ correct on a simple open-set word recognition test (Multisyllabic Lexical Neighborhood Test) administered using monitored live voice (70 dB SPL). In older children (≥ 4 years of age), lack of hearing aid benefit is defined as scoring $\leq 12\%$ on a difficult open-set word recognition test (Phonetically Balanced-Kindergarten Test) or $\leq 30\%$ on an open-set sentence test (Hearing In Noise Test for Children) administered using recorded materials in the sound field (70 dB SPL).*

Intended Purpose

The Naída CI M90 sound processor is an accessory of an auditory active implantable system, the HiResolution Bionic Ear system. The HiResolution Bionic Ear system is intended to provide auditory sensation via electrical stimulation of the auditory nerve for individuals with severe to profound bilateral or unilateral sensorineural hearing loss. Severe hearing loss is defined as audiometric thresholds greater than or equal to 70 dB HL, but less than 90 dB HL. Profound hearing loss

is defined as audiometric thresholds greater than or equal to 90 dB HL. The Naída CI M90 sound processor is a behind the ear (BTE) sound processor which works together with the implant to bypass the damaged part of the inner ear and converts sound picked up by the microphone or streamed via wireless communication into electrical signals that are used by the cochlear implant to enable hearing. The Naída CI M90 is the premium version with full access to multiple automatic programs and features, including bimodal and bilateral and supports compatibility with acoustic amplification.

Intended Users

The intended users of the Naída CI M90 sound processor are recipients of Advanced Bionics cochlear implants, their caregivers if applicable, and hearing care professionals.

The recipient, or their caregiver, should, at a minimum, be capable of changing the battery, connecting the earhook, placing and removing the sound processor from their ear, and placing and removing the headpiece from the implant site.

The hearing care professional should be trained in the use and fitting of the cochlear implant system.

Intended Use Environment

The intended use environments of the Naída CI M90 sound processor are healthcare and daily living environments. The Naída CI M90 sound processor is designed to be worn daily during the user's normal waking hours.

Compatibility

The Naída CI M90 is compatible with the following implant types from Advanced Bionics:

- **CII**
- **HiRes™ 90K**
- **HiRes™ 90K Advantage**
- **HiRes™ Ultra**
- **HiRes™ Ultra 3D**

The Naída CI M90 is compatible with the following Advanced Bionics products:

- **M Batteries, M Zn-Air Battery Pak¹, M Waterproof Battery**
- **M T-Mic™ microphone, M Earhook, M Acoustic Earhook**
- **Slim HP, Slim HP Mic, Slim HP AquaMic™ headpiece**
- **M Wearing Clip, M Snuggie™ accessory, M Retention Cuff, M Listening Check™ accessory**

The Naída CI M90 is compatible with the following programming interfaces:

- **Advanced Bionics Clinical Programming Interface CPI-3 with the M Programming Cable**
- **NoahLink Wireless programmer**

The Naída CI M90 is compatible with the following accessory products:

- **AB Remote mobile application**
- **Phonak Roger transmitters**

- **Phonak RemoteControl**
- **Phonak TV Connector**
- **Phonak PartnerMic**
- **Phonak Naída™ Link M hearing instrument**

The Naída CI M90 sound processor is Bluetooth 4.2 certified. Any mobile device that supports Bluetooth, Hands Free Profile (HFP) and/or Advanced Audio Distribution Profile (A2DP) should work with your Naída CI M90 sound processor. Bluetooth compatibility information should be used as a reference guide only, because hardware and software levels and releases may vary, and numerous hardware and software combinations are possible. Advanced Bionics assumes no responsibility for the compatibility; thus, it is the user's responsibility to do product compatibility testing by himself before making purchasing decisions on smartphones and other Bluetooth-enabled devices.

Operating Instructions

Powering Naída CI M90

The sound processor is turned on when a charged battery is attached to the sound processor. When the battery is engaged, the orange LED located in the middle of the multifunction button will indicate battery charge.

The startup program is selected by the cochlear implant professional during programming.

To power down the Naída CI M90 sound processor, simply remove the battery cartridge.

Attaching the Battery

Hold the Naída CI M90 in one hand. Rotate the headpiece cable upward in order to avoid bumping it upon battery placement.

Align the battery connector with the sound processor housing side of the connector.

Slide the battery onto the sound processor until it clicks into place. Do not force the battery onto the sound processor. The batteries are designed to be inserted in only one direction; applying force may damage the equipment.

Removing the Battery

Hold the Naída CI M90 in one hand. In the other hand, hold the battery. Rotate the headpiece cable upward in order to avoid bumping it upon battery removal.

Firmly slide the battery away from the M T-Mic, M Earhook, or M Acoustic Earhook.

Continue sliding the battery until it separates from the sound processor.

Approved Naída CI M90 Power Supplies

POWER SUPPLY	NOMINAL VOLTAGE	TYPE	RATED CAPACITY
M Battery	3.6V – 3.7V	Li-ion (rechargeable)	0.3Wh – 0.8Wh
M Waterproof Battery	3.7V	Li-ion (rechargeable)	0.8Wh

POWER SUPPLY	NOMINAL VOLTAGE	TYPE	RATED CAPACITY
M Zn-Air Battery Pak ¹	2.3V	Zn-air (disposable)	1.4Wh

Connecting the earhook to the sound processor

There are three earhook options for the Naída CI M90 sound processor: the M T-Mic, the M Earhook, and the M Acoustic Earhook.

Beginning with the earhook disconnected, align the earhook flush with the sound processor. Hold the earhook firmly at the top and push towards the sound processor until it meets flush with the sound processor.

Using the provided pin removal tool, place one of the pins in the small hole on the side of the tool. This will help to steady the pin as you line it up with the pin hole on the earhook.

Line up the pin with the hole in the earhook and gently push it through the hole. The pin will move from one side of the sound processor, through the earhook to the other side of the sound processor, holding the earhook in place. The pin should not come out the other side.

If a small portion of the pin extends beyond the sound processor, use the side of the tool to push the pin inwards so it lies flush with the sound processor.

Pull the tool away from the sound processor, leaving the pin inside the sound processor.

Removing the earhook

Align the pointed end of the tool with the hole on the sound processor that contains the pin holding the earhook on the sound processor.

Gently press the tool in the hole, pushing the pin out the other side.

The pin does not need to be removed from the sound processor. As long as the tool is pushed in flush on one side of the sound processor, the pin should be displaced enough to adequately remove the earhook. Remove the tool from the sound processor. Gently grasp the earhook at the top near where it connects with the sound processor and pull forward until it is separated from the sound processor.

Attaching the headpiece cable to the sound processor

To attach the headpiece cable to the sound processor, hold the headpiece cable by the cable's strain relief. Align the connector on the cable with the connector on the sound processor and gently push the cable into the sound processor until it snaps into place.

Removing the Slim HP cable from the sound processor

Remove the battery from the processor before removing the cable from the processor. To remove the cable from the sound processor, always hold the cable's strain relief and gently pull it away.

NOTE: To maximize cable life, the cable should only be unplugged from the sound processor when changing wearing configuration or when replacing the cable.

Naída CI M90 LEDs

The LED is a programmable feature that provides visual information about the Naída CI M sound processor status, battery life, program position, and sound processor error conditions.

INDICATION	Battery Status at Start-Up
COLOR	Orange
BEHAVIOR	Blinks at start-up with use of rechargeable batteries (not available with M Zn-Air Battery Pak¹): • 4 quick blinks indicate that the battery is fully-charged. • 2 or 3 quick blinks indicate that the battery is sufficiently charged to power the sound processor. • 1 quick blink indicates that the battery is nearly depleted. • No blinking indicates the battery is depleted; replace with a charged or new battery.
PROGRAMMABLE	No

INDICATION	Low Battery
COLOR	Orange
BEHAVIOR	Continuous Blinking
PROGRAMMABLE	No

INDICATION	Current Program
COLOR	Green
BEHAVIOR	Blinks at start-up after battery status and upon program change • 1 long blink is AutoSense program • 1 short blink is program 1 • 2 short blinks is program 2 • 3 short blinks is program 3 • 4 short blinks is program 4

PROGRAMMABLE	Yes
---------------------	-----

INDICATION	Flight Mode
COLOR	Orange
BEHAVIOR	One long blink at start-up after battery and program status
PROGRAMMABLE	No

INDICATION	Progressive Level Change
COLOR	Purple
BEHAVIOR	1 blink upon changing
PROGRAMMABLE	Yes

INDICATION	Loud Input
COLOR	Green
BEHAVIOR	Blinks during use in response to loud input sounds. Indicates the sound processor and microphone are responding to sound.
PROGRAMMABLE	Yes

INDICATION	Streaming Input
COLOR	Blue
BEHAVIOR	Blinks while streaming audio to the sound processor. Indicates the sound processor and microphone are streaming audio.
PROGRAMMABLE	Yes

INDICATION	Loss of Lock with the Implant
-------------------	-------------------------------

COLOR	Red
BEHAVIOR	Slow blinking (once per second)
PROGRAMMABLE	Yes

INDICATION	Wrong Implant
COLOR	Red
BEHAVIOR	Fast blinking (more than once per second)
PROGRAMMABLE	No

INDICATION	Error Condition (remove and reattach battery to reset sound processor)
COLOR	Red
BEHAVIOR	Solid Red
PROGRAMMABLE	No

INDICATION	Detect Sound Processor
COLOR	Green
BEHAVIOR	3 blinks upon command from mobile application
PROGRAMMABLE	No

Naída CI M90 Internal Alarms

The internal alarms are programmable and provide auditory information about the Naída CI M90 sound processor. Your cochlear implant professional can adjust the volume and pitch of the internal alarms to suit your preference.

INDICATION	BEHAVIOR	PROGRAMMABLE
Low Battery	2 beeps	No
Program Change	Jingle is AutoSense program • 1 beep is program 1 • 2 beeps is program 2 • 3 beeps is program 3 • 4 beeps is program 4	Yes
Volume Change	• Long beep at end of range • 2 short beeps at middle of range • 1 short beep at intermediate steps	Yes
Pairing Successful	Ascending beeps	Yes
TV Connector Available	Ascending beeps	Yes
Phone Ring Tones	Ringing beeps	No

Using the Naída CI M90 Multi function Button

The multifunction button has several functions which can be programmed by the hearing care professional. The button can function as a volume control and/or a program change depending on the programming. Please ask your cochlear implant professional to confirm how your sound processor is programmed.

SITUATION	SHORT PRESS	LONG PRESS
Standard Use	Volume up or down	Program change
During Audio Streaming	Streaming volume up or down	Program change
During Phone Calls	Accept phone call (while ringing) or phone call volume up or down (during call)	Reject phone call (while ringing) or end phone call (during call)

Connectivity Overview

The Naída CI M90 sound processor can be connected to Bluetooth-enabled devices for phone calls, VoIP calls, audio streaming, or for use with the AB Remote mobile application.

Pairing and Connecting the Naída CI M90 to a Bluetooth-Enabled Device

- 1. Ensure that Bluetooth wireless technology is enabled on your device (e.g., phone, tablet) and search for Bluetooth-enabled devices in the connectivity settings menu.**
- 2. Power on the Naída CI M90 sound processor by connecting the sound processor to a power source. The sound processor will be discoverable by your Bluetooth-enabled device for up to three minutes, or until it is successfully paired to a device.**

- 3. Your device should show a list of Bluetooth-enabled devices. Select the sound processor from the list to pair your sound processor. If you use two Naída CI M90 sound processors, selecting one sound processor for pairing will simultaneously pair both sound processors. A beep confirms successful pairing.**
- 4. The connection will be maintained as long as the device remains on and within range.**

After your sound processor has been paired to your Bluetooth-enabled device, the sound processor will automatically connect again when switched on. Your sound processor can be paired to two Bluetooth-enabled devices at a time. The paired sound processor can be unpaired from the Bluetooth device by utilizing the Bluetooth management menu of your device.

NOTE: Once a sound processor has been paired to a device it is recommended to keep paired devices in silent mode to avoid dropouts related to device notifications.

Managing Bluetooth Phone Calls with the Naída CI M90

When the sound processor is directly connected with a Bluetooth-enabled phone, you will be able to hear the caller's voice directly through the sound processor. The sound processor picks up your voice through the sound processor microphones.

When making a phone call, you will hear a dialing tone through the sound processor. The sound processor microphones pick up your voice. When receiving a call, a call notification will be heard through the sound processor. You can accept the call with a short press on the upper or lower part of the multifunction button on the sound processor or directly on the phone.

You can reject or end a call with a long press on the upper or lower part of the multifunction button on the sound processor or directly through the phone.

Streaming Audio via Bluetooth with the Naída CI M90

After successful pairing of your processor(s) to your Bluetooth-enabled device, ensure that your device is not muted and volume is not set to minimum. Then, simply start the audio on the Bluetooth-enabled device and the sound processor will start streaming the audio. Environmental balance and access to surrounding sounds can be managed using the multifunction button, the AB Remote mobile application, or the Phonak RemoteControl.

Pairing the Naída CI M90 to Phonak Accessories

Your Naída CI M90 can be paired to the following Phonak accessories:

- Phonak RemoteControl
- Phonak TV Connector
- Phonak PartnerMic
- Phonak Roger transmitters
- Phonak Naída Link M hearing instrument

Please refer to the accessory user guide for initial setup and pairing instructions.

Switching between Multiple Audio Sources with the Naída CI M90

A phone call from a Bluetooth-enabled device has highest priority and will pause other audio streaming from a connected device. In order to switch audio streaming source, simply pause the audio from the current source and start streaming from another source.

Placing the Naída CI M90 in Flight Mode

The Naída CI M90 communicates wirelessly with other devices in the 2.40 GHz to 2.48 GHz frequency range. When flying, some operators require all devices to be switched into flight mode. Entering flight mode will not disable normal sound processor functionality but will disable Bluetooth connectivity functions.

Activating and Deactivating Flight Mode

1. Hold down either the top or bottom of the multifunction button while attaching the battery.
2. Continue to hold down the button while the sound processor starts up, until a long, solid orange LED is visible, approximately 10 seconds. Note that the start-up LEDs that indicate the battery life of the rechargeable battery and the program number will be seen before the long, solid orange flight mode LED.
3. Removing and reattaching the battery will take the sound processor out of flight mode.

Care and Maintenance

Storing the Naída CI M90

- When not in use, store the Naída CI M90 in the equipment case provided.
- If product is stored at the minimum or maximum storage temperature, bring the product into an ambient temperature environment (approximately 20°C (68°F)) for 3 minutes before using the product.

Recommended Operating and Storage Temperature, Humidity and Pressure Ranges

CONDITION	MINIMUM	MAXIMUM
Operating Temperature	0°C (32°F)	45°C (113°F)*
Storage Temperature	-20°C (-4°F)	55°C (131°F)
Relative Humidity	0%	95%
Pressure	70 kPa	106 kPa

*This product may reach a temperature of 49°C when used in a 45°C environment.

Cleaning and Maintenance

Wipe the outside of the Naída CI M90 with a soft, dry cloth. Do not immerse in liquids or use cleaning agents for cleaning the sound processor. If the Naída CI M90 sound processor is exposed to moisture, place the device in a dryer before next use. Before using hair spray or applying cosmetics you should remove the sound processor from the ear because these products may cause damage.

Expected Life

The expected life of the Naída CI M90 sound processor is 5 years.

Ingress Protection (IP) Ratings

The following product configuration has an IP rating of 22 (protection against insertion of solid objects ≥ 12.5 mm diameter; protection against failure due to dripping water when tilted up to 15°):

- **The Naída CI M90 sound processor with the M Listening Check, an M Battery or M Zn-Air Battery Pak¹, a M T-Mic, M Earhook, or M Acoustic Earhook, and a headpiece.**

The following product configuration has an IP rating of 52 (protection against dust; protection against failure due to dripping water when tilted up to 15°):

- **The Naída CI M90 sound processor inside the M Waterproof Battery with a non-waterproof headpiece.**

The following product configurations have an IP rating of 54 (protection against dust; protection against failure due to water splashes from all directions):

- **The Naída CI M90 sound processor with an M Battery or M Zn-Air Battery Pak¹, an M T-Mic or M Acoustic Earhook, and a Slim HP or Slim HP Mic.**
- **The Naída CI M90 sound processor with an M Zn-Air Battery Pak¹, an M Earhook, and a Slim HP or Slim HP Mic.**

The following product configuration has an IP rating of 57 (protection against dust; protection against failure due to one-time immersion for 30 minutes to a depth up to 1 meter and after drying overnight in a drying system):

- **The Naída CI M90 sound processor with an M Battery, an M Earhook, and a Slim HP or Slim HP Mic.**

The following product configuration has an IP rating of 68 (complete protection against dust penetration; protection against failure due to continuous immersion in water up to 3 meters):

- **The Naída CI M90 sound processor inside the M Waterproof Battery with Slim HP AquaMic.**

DISPOSAL AND HAZARDOUS MATERIALS

Dispose of the Naída CI M90 sound processor and accessories in accordance with applicable national and local regulations. Hazardous materials: not applicable to this product.

AVAILABLE PRODUCT OPTIONS AND MODEL NUMBERS

PRODUCT OPTION	MODEL NUMBER
Sand Beige	CI-5293-120
Chestnut	CI-5293-130
Silver Grey	CI-5293-140

PRODUCT OPTION	MODEL NUMBER
Velvet Black	CI-5293-150
Alpine White	CI-5293-110
Resin Beige	CI-5293-240

CLINICAL BENEFITS

The Intended Clinical Benefit of the Naída CI M90 sound processor as part of the HiResolution Bionic Ear System is to:

- *provide useful hearing to individuals with severe-to-profound hearing loss via electrical stimulation of the auditory nerve.*
- *provide a combined type of stimulation, electric plus acoustic, when sound processors are enabled with the Acoustic Earhook.*

Advanced Bionics has performed two clinical studies of AutoSense OS on the Naída CI M90 sound processor:

- *a confirmatory study comparing a Naída CI M sound processor to a Naída CI Q sound processor with 10 subjects using T-Mics and electrical stimulation only.*
- *a supplemental study comparing AutoSense OS On with AutoSense OS Off with 10 subjects using a Naída CI M sound processor with either the M Acoustic Earhook or the M T-Mic.*

Confirmatory Study

A total of 10 subjects were enrolled in a prospective, confirmatory study using a within-subjects repeated-measures design where each subject served as his/her own control. All 10 subjects used the used the M T-Mic and electrical stimulation only. The data presented below show the comparison between AutoSense OS on a Naída CI M sound processor to AutoSound on a Naída CI Q sound processor in quiet and noise as well as a comparison of AutoSense OS On vs AutoSense OS Off on a Naída CI M sound processor in noise.

Sentences in Quiet at 65 dBA

Similar sentence recognition scores in quiet were observed between AutoSense OS on a Naída CI M sound processor and AutoSound on a Naída CI Q sound processor. See the table below.

	NAÍDA CI Q 2 LIST AUTO SOUND AVERAGE SCORE (%)	NAÍDA CI M 2 LIST AUTO SENSE OS AVERAGE SCORE (%)
N	10	10
Mean (SD)	88.31 (5.817)	87.30 (10.214)
Median	89.05	90.70
Min, Max	80.6 97.4	72.4 98.6

Sentences in Noise at 65 dBA, + 5 dB SNR Multi-Talker Babble

Sentence recognition scores in noise with AutoSense OS on a Naída CI M sound processor were better than sentence scores with AutoSound on a Naída CI Q sound processor and also were better than sentence scores with the AutoSense OS Off program on a Naída CI M sound processor. See table below.

	NAÍDA CI Q	NAÍDA CI M	
	2 LIST AUTO SOUND AVERAGE SCORE (%)	2 LIST AUTO SENSE OS AVERAGE SCORE (%)	2 LIST AUTO SENSE OS OFF AVERAGE SCORE (%)
N	10	10	10
Mean (SD)	52.95 (31.163)	76.37 (19.078)	42.93 (30.801)
Median	48.25	79.50	39.15
Min, Max	13.4 94.7	31.8 92.8	9.3 91.2

Statistical analysis of the sentence recognition data shows that AutoSense OS on a Naída CI M sound processor performs better than AutoSound on a Naída CI Q sound processor or the AutoSense OS Off program on a Naída CI M sound processor when tested in noise. When tested in quiet, statistical analysis of the data shows that AutoSense OS on a Naída CI M sound processor is no worse than AutoSound on a Naída CI Q sound processor.

Self-Assessment Questionnaire Results

Overall subjects' responses showed higher ratings after using AutoSense OS on a Naída CI M sound processor for 2-3 weeks as compared to ratings on their own processor at the beginning of the study, indicating better acceptability of AutoSense on a Naída CI M processor. All subjects agreed or strongly agreed that AutoSense on a M processor is acceptable for use.

Supplemental Study

A total of 10 subjects were enrolled in a prospective, supplemental study using a within-subjects repeated-measures design where each subject served as his/her own control. Five subjects used the M Acoustic Earhook, while the remaining 5 used the M T-Mic and electrical stimulation only. The data presented below are pooled between the two cohorts.

Sentences in Quiet at 65 dBA

Across all subjects, similar sentence recognition scores in quiet were observed between AutoSense OS and the AutoSense OS Off program. See the table below.

	2 LIST AUTOSENSE OS AVERAGE SCORE (%)	2 LIST AUTOSENSE OS OFF AVERAGE SCORE (%)
N	10	10
Mean (SD)	88.92 (11.068)	89.76 (11.677)
Median	94.65	92.88

	2 LIST AUTOSENSE OS AVERAGE SCORE (%)	2 LIST AUTOSENSE OS OFF AVERAGE SCORE (%)
Min, Max	63.2, 96.7	61.2, 99.3

Sentences in Noise at 65 dBA, + 5 dB SNR Multi-Talker Babble

Across all subjects, better sentence recognition scores in noise were observed when using AutoSense OS as compared to the AutoSense OS Off program. See table below.

	2 LIST AUTOSENSE OS AVERAGE SCORE (%)	2 LIST AUTOSENSE OS OFF AVERAGE SCORE (%)
N	10	10
Mean (SD)	79.29 (16.439)	54.19 (23.826)
Median	86.98	57.53
Min, Max	45.8, 95.3	19.8, 81.4

In quiet, 100% of subjects' AutoSense OS scores were within 10% of their AutoSense OS Off scores and in noise 90% of subjects' AutoSense OS scores were better than their AutoSense OS Off scores by 10% or more. These observed trends further demonstrate that automatic activation of the AutoSense OS features does not negatively impact sentence recognition in quiet or noise conditions. Additionally, the improvement in sentence recognition in multi-talker babble noise with AutoSense OS is evident.

Follow-Up Self-Assessment Questionnaire Results

Results were similar across sound quality and listening comfort ratings. All 10 subjects (100.0%) reported acceptable sound quality in quiet, with 9 subjects (90.0%) reporting strong agreement and 1 subject (10.0%) reporting slight agreement. Nine subjects (90.0%) reported acceptable sound quality in noise (5 subjects reported "Strongly Agree," 4 subjects reported "Slightly Agree"), with 1 subject reporting slight disagreement. All 10 subjects (100.0%) strongly agreed that listening in quiet was comfortable. Nine subjects (90.0%) agreed that listening in noise was comfortable (7 subjects reported "Strongly Agree," 2 subjects reported "Slightly Agree"), with 1 subject reporting slight disagreement.

All 10 subjects (100.0%) reported that AutoSense OS on the sound processor was acceptable for use and daily wear, with 8 subjects (80.0%) reporting strong agreement and 2 subjects (20.0%) reporting slight agreement with the statements "The sound processor is acceptable for use" and "The sound processor is acceptable for daily wear". The majority of subjects (80.0%) also agreed that AutoSense OS on the sound processor met their listening needs. The 2 subjects that slightly disagreed reported challenges in noisy environments, although one of these subjects agreed that sound quality was acceptable and listening was comfortable in noise. These questionnaire results indicate that the sound processor is acceptable for use and satisfies users' listening needs.

In the United States, ClearVoice is approved for pediatric use in children 6 years and above who 1) are able to complete objective speech perception testing in order to determine speech performance and 2) are able to report a preference for different coding strategies or features.

ClearVoice is only available in markets where ClearVoice has received regulatory approval. Contact Advanced Bionics for more information.

COUNSELLING AND FITTING RECOMMENDATIONS

- **WindBlock:** The intent of the WindBlock sound processing feature is to reduce wind noise arriving at the microphone to improve comfort and ease of listening when listening in the presence of wind. This is an automatic feature that is enabled in the AutoSense program and can be enabled in other programs.
- **SoundRelax:** The intent of the SoundRelax sound processing feature is to quickly respond to sudden and/or unexpected sounds. This is an automatic feature that is enabled in the AutoSense program and can be enabled in other programs.
- **EchoBlock:** The intent of the EchoBlock sound processing feature is to improve comfort and ease of listening when listening in environments with a high degree of reverberation or echo. This is an automatic feature that is enabled in the AutoSense program and can be enabled in other programs.

- **UltraZoom:** The intent of the UltraZoom sound processing feature is to improve the signal to enhance listening when focused on a person speaking in front of you. In the AutoSense program, the feature will automatically enter or exit this mode depending on the surrounding noise in the environment. UltraZoom can be enabled in a static mode, which is referred to as Fixed Directional mode.
- **StereoZoom:** The StereoZoom feature is designed only for use with bilateral recipients. The intent of this feature is to use a bilateral beamformer to create a narrow beam to focus on one person in front of you when in a noisy environment. This is an automatic feature that is enabled in the AutoSense program and can be enabled in other programs.
- **RealEar Sound:** The intent of the RealEar Sound sound processing feature is to simulate the natural sound quality of the M T-Mic microphone. RealEar Sound is active in quieter listening environments when the M T-Mic microphone is not in use.
- **WhistleBlock:** The intent of the WhistleBlock sound processing feature is to minimize acoustic feedback that can occur when the M Acoustic Earhook is used with the Naída CI M90 sound processor. This is an automatic feature.
- **NoiseBlock:** The intent of the NoiseBlock sound processing feature is to improve comfort when listening in noisy environments. This is an automatic

feature that is available when the M Acoustic Earhook is in use. It is enabled in the AutoSense program and can be enabled in other programs.

- **AutoSense OS:** The intent of the AutoSense program is to classify the listening environment and activate sound processing features based on the listening environment. The AutoSense program is enabled by default as the start-up program.

Cautions

- UltraZoom and StereoZoom are beamforming algorithms; therefore, by design they focus to the front and dampen sound from the side and back of the recipient.
- The Naída CI sound processor is designed so it can be worn in an off-ear wearing configuration. We do not recommend the processor be worn in the off-ear wearing configuration while using any beamforming algorithm.
- The intent of the WindBlock, EchoBlock, SoundRelax, and NoiseBlock algorithms is to attenuate a portion of the signal in order to provide comfort in specific environments. Attenuation of the signal may affect sound quality.

Note: In the United States, Fixed Directional mode, StereoZoom, UltraZoom, Sound Relax, NoiseBlock, WindBlock, and EchoBlock are approved for use in pediatric recipients 6 years and above who are 1) able to complete objective speech perception testing in order to determine speech performance and 2) able to report a preference for different coding strategies or features.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

*A summary of safety and clinical performance (SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI:- 08400944CI5293YE
<https://ec.europa.eu/tools/eudamed>*

Users of the Naída CI M90 sound processor in the European Union should report any serious incident to their local competent authority and to Advanced Bionics.


GUIDELINES FOR IT SECURITY

Intended Use Environment

The Naída CI M90 sound processor system is intended for use in healthcare and home environments. The home environment is extended to include use outdoors and during travel (e.g., airplanes). The Naída CI M90 sound processor system is intended to operate with the AB Remote application, the Phonak RemoteControl, the

Target CI fitting software, and with any Bluetooth-enabled device used for audio streaming. The user is responsible for security of connections with other devices, including pairing the sound processor with other devices in a secure environment.

System and Network Diagram

System Interfaces

COMPONENT	Sound Processor
INTERFACE	Bluetooth
DATA TRANSFER DIRECTION	Bidirectional

COMMUNICATION PROTOCOL	Bluetooth Classic: Hands-Free Profile with wideband speech option v1.6 (SPP v1.2, RFCOMM v1.2 over L2CAP) Bluetooth Classic: Advanced Audio Distribution Profile v1.3 (AVDTP v1.3 over L2CAP) Bluetooth Classic: Audio/Video Remote Control Profile v1.5 (AVCTP v1.4 over L2CAP)
SECURITY	Bluetooth Classic: E0 encryption Bluetooth LE: AES-CCM 128 encryption

COMPONENT	Sound Processor
INTERFACE	RF
DATA TRANSFER DIRECTION	Bidirectional
COMMUNICATION PROTOCOL	Advanced Bionics proprietary
SECURITY	None

COMPONENT	Remote Control
INTERFACE	Bluetooth
DATA TRANSFER DIRECTION	Bidirectional
COMMUNICATION PROTOCOL	Bluetooth LE: Proprietary GATT (Generic Attribute Profile) over L2CAP
SECURITY	Bluetooth LE: AES-CCM 128 encryption

COMPONENT	Roger Devices
------------------	---------------

INTERFACE	Roger DMI
DATA TRANSFER DIRECTION	Bidirectional
COMMUNICATION PROTOCOL	Phonak proprietary
SECURITY	32-bit code defining address and hopping sequence

COMPONENT	TV Connector PartnerMic
INTERFACE	Airstream
DATA TRANSFER DIRECTION	Bidirectional
COMMUNICATION PROTOCOL	Phonak proprietary
SECURITY	AES-CTR 128 encryption

Cybersecurity Configuration

There is no specific Naída CI M90 sound processor configuration needed to help ensure cybersecurity. Ensure pairing with Bluetooth devices is completed in a secure environment.

Integrated Cybersecurity Functions

The Naída CI M90 system does not include user authorization or authentication mechanisms. The sound processor is configured to connect with a single Advanced Bionics cochlear implant when fit with the Target CI fitting software. All events are logged to an internal event log. Advanced Bionics can interpret and analyze the log if the

device is returned. Data integrity and quality of service is ensured through sound processor and implant error correction and detection mechanisms.

System, Configuration, and Data Backup/Restore

The user should visit his or her cochlear implant professional if any unintended behavior of the sound processor is observed.

Cybersecurity Incident Response

The Naída CI M90 does not detect cybersecurity attacks or intrusions. In case of abnormal behavior of the sound processor, please visit your cochlear implant professional or contact Advanced Bionics.

Patches and Updates

Advanced Bionics does not patch the Naída CI M90 sound processor. Firmware updates may be made available from time to time and will require visiting your cochlear implant professional to reimage your sound processor. Authenticity of firmware can be verified by your cochlear implant professional using the Target CI fitting software.

Available Training

Users of the Naída CI M90 sound processor are trained by their cochlear implant professional at the initial fitting session.

End of Cybersecurity Support

Advanced Bionics products are supported until they no longer have regulatory approval in the country of use or until the product is obsolesced by Advanced Bionics.

GUIDANCE AND MANUFACTURER'S DECLARATION

Per IEC 60601-1-2

Electromagnetic Emissions

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The Naída CI M90 sound processor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Naída CI M90 sound processor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Electromagnetic Immunity

The Naída CI M90 sound processor is intended for use in the electromagnetic environment specified below. The customer or the user of the Naída CI M90 sound processor should assure that it is used in such an environment.

IMMUNITY	IEC 60601 TEST LEVEL	COM-PLIANCE LEVEL ^a	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electro-static discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. As with the operation of other electronic devices, precaution should be taken to not generate ESD.
Power frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

IMMUNITY	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL ^a	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz 3	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the Naída CI, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2\sqrt{P} < 800 \text{ MHz}$ $d = 2.3\sqrt{P} \geq 800 \text{ MHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^b should be less than the compliance level in each frequency range.^c Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

- a. Essential performance of the Naída CI per IEC 60601 requirements is defined as auditory stimulation within safe amplitudes.
- b. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Naída CI is used exceeds the applicable RF compliance level above, the Naída CI should be observed to verify normal operation.
- c. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Separation Distances between RF Communications Equipment and Naída CI M90

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE NAÍDA CI M90

The Naída CI is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Naída CI can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Naída CI as recommended below, according to the maximum output power of the communications equipment.

RATED MAXIMUM OUTPUT POWER OF TRANSMITTER (W)	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER (M)	
	$d = 1.2\sqrt{P}$ < 800 MHz	$d = 2.3\sqrt{P}$ ≥ 800 MHz
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		

Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Quality of Service

The proprietary wireless communication protocol used to connect to Phonak wireless accessories meets a Quality of Service of an audio error rate of 0.008 which translates to a maximum of 2 missed audio frames per second and maintains good audio quality.

The Sonova Binaural Protocol used to communicate between hearing instruments meets a Quality of Service of a frame error rate (FER) of 1.6% to maintain good audio quality.

All Bluetooth-related features of the product such as HFP, A2DP, and BLE services meet Quality of Service requirements that are determined by the Bluetooth Special Interest Group (SIG).

Wireless Safety

Wireless Safety Measures have been put into place for the Sonova proprietary protocols:

- A proprietary network protocol which creates and manage the devices in the network.
- A unique binaural group ID programmed to the processor and its accessories at fitting time for SBP which is validated during each transmission.
- A unique pairing key established at pairing time for DM devices which is validated during transmission.
- A cyclic redundancy check (CRC) which ensures the integrity of remote-control and audio data.
- A system architecture that isolates the processor's critical functionality of stimulation parameters from the wireless functionality.

Bluetooth safety features have been implemented according to the Bluetooth SIG requirements and tested according to those specifications.

Radio Information for your Sound Processor

Antenna type	Magnetic loop antenna
Operation frequency	2.4 GHz – 2.48 GHz
Modulation	GFSK
Radiated power	<2.5mW
Bluetooth	
Range	~1m
Bluetooth	4.2 Dual Mode
Profiles supported	HFP (Hands Free Profile), A2DP

This instrument is certified under:

- **FCC ID: 2AU60-ABBTE2**
- **IC: 25853-ABBTE2**

This device complies with Part 15 of the FCC Rules and with RSS-247 of Industry Canada. Operation is subject to the following two conditions:

- 1. This device may not cause harmful interference, and**
- 2. This device must accept any interference received, including interference that may cause undesired operation.**

Changes or modifications made to this equipment not expressly approved by Advanced Bionics may void the FCC authorization to operate this equipment. This Class B digital apparatus complies with Canadian ICES-003.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- **Reorient or relocate the receiving antenna.**
- **Increase the separation between the equipment and receiver.**
- **Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.**
- **Consult your cochlear implant professional or an experienced radio/TV technician for help.**

†Note: Certain products are only available in markets that have received regulatory approval. Contact Advanced Bionics for more information.

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Advanced Bionics LLC
Advanced Bionics GmbH