ADVANCED
BIONICS

POWERFUL CONNECTIONS

Slim HP, Slim HP Mic, Slim HP AquaMic™

CI-5320, CI-5322, CI-5321

Slim HP Color Caps

CI-7131, CI-7132

instructions for use

English, Français, Español, Deutsch, Nederlands, Italiano

Contents

PURPOSE AND SCOPE OF IFU	4
LIMITATIONS AND CONTRAINDICATIONS	5
PRECAUTIONS, CAUTIONS, AND WARNINGS	5
UNDESIRABLE SIDE EFFECTS	7
PRODUCT DESCRIPTION	7
INDICATIONS FOR USE	10
INTENDED PURPOSE	11
<i>Intended Users</i>	13
<i>Intended Use Environment</i>	14
COMPATIBILITY	14
OPERATING INSTRUCTIONS	14
<i>Assembly Instructions</i>	14
<i>Cables</i>	15
<i>Color Caps</i>	17
<i>Magnets</i>	18
<i>For use of Slim HP AquaMic in water</i>	18
CARE AND MAINTENANCE	18
<i>Storing the Slim HPs</i>	18
<i>Cleaning and Maintenance</i>	18
<i>Recommended Operating and Storage Temperature Ranges</i>	19
<i>Expected Life</i>	19
<i>Ingress Protection (IP) Ratings</i>	19

DISPOSAL AND HAZARDOUS MATERIALS	20
AVAILABLE PRODUCT OPTIONS AND MODEL NUMBERS	20
<i>Slim HP</i>	20
<i>Slim HP Mic</i>	21
<i>Slim HP AquaMic</i>	21
<i>Slim HP Color Caps</i>	22
CLINICAL BENEFITS	24
SUMMARY OF SAFETY AND CLINICAL PERFORMANCE	24

SYMBOLS TABLE

Labeling symbols and their meanings

	European Community Mark of Conformity. Authorized to affix the CE Mark in 2021
	Date of Manufacture
	Manufacturer
	Model number
	LOT number
	Medical Device
	Unique Device Identifier
	European authorized representative
	Type of Protection: BF
	Suitable temperature range for transport and storage
	See Instructions for Use

	Fragile
	Do not get wet
	Suitable relative humidity range for transport and storage
	Caution (refer to “Cautions & Warnings” section for further information)
	Dispose of in accordance with applicable national and local regulations
	MR Unsafe
IPxx	Ingress Protection Rating

PURPOSE AND SCOPE OF IFU

This IFU is designed to help cochlear implant recipients, and their caregivers if applicable, to understand the use and care of the Slim HP, Slim HP Mic, Slim HP AquaMic, and Slim HP Color Caps. It is also meant to be used by hearing care professionals trained in the fitting of cochlear implants to counsel the recipients on the use of the Slim HP, Slim HP Mic, Slim HP AquaMic, and Slim HP Color Caps.

Limitations and Contraindications

As the Slim HP, Slim HP Mic, Slim HP AquaMic, and Slim HP Color Caps are components of the HiResolution™ Bionic Ear System, the following contraindications stated for the HiResolution Bionic Ear System are applicable: Deafness due to lesions of the acoustic nerve or central auditory pathway; cochlear ossification that prevents electrode insertion; absence of cochlear development.

Precautions, Cautions, and Warnings

- *This device should be used only by the individual for whom it is prescribed.*
- *Do not remove or replace the Slim HP cable from the sound processor or waterproof wearing solution while underwater. If the cable becomes accidentally detached during water use, remove the headpiece and sound processor from the water and dry both devices prior to re-connecting the cable to the sound processor.*
- **CHOKING HAZARD:** *Contains small parts that pose a hazard of inhalation, choking, or ingestion. If any parts are swallowed or inhaled, consult a physician or hospital immediately.*
- *Do not leave children unattended with or allow children to play with the Slim HPs or Slim HP Magnets.*
- *Do not allow children to place the Slim HPs or Slim HP Magnets in the mouth.*
- *The large Slim HP Color Cap is for use by recipients older than 3 years of age, as this color cap and the Slim HP magnet may become displaced upon dropping.*

- *Due to increased risk of strangulation, do not use cables longer than 81cm (32") with young children.*
- *Remove your sound processor and consult your cochlear implant professional if uncomfortable sounds are heard or in case of discomfort, pain or skin irritation.*
- *Be sure to choose a cable length that does not require looping or bundling of the cable due to excess length.*
- *The headpiece must only be used by the intended person and not by others.*
- *Do not use a headpiece, magnet, or any associated part, that is broken, damaged, or cracked.*
- *The color cap must be attached when the headpiece is in use.*
- *Store spare magnets safely and away from cards that may have a magnetic strip (e.g. credit cards) and from other spare magnets or headpieces.*
- *High levels of static electricity may damage electronic components of your sound processor or implant. Care should be taken to avoid exposure of the system to static electricity.*
- *Use the lowest strength magnet required to ensure sufficient magnetic attraction of the headpiece.*
- *Contact your cochlear implant professional to adjust the magnet strength if the magnetic attraction of the headpiece is insufficient.*
- *If the headpiece magnet is too strong, pressure sores may develop at the headpiece site. If this happens,*

or if you experience any discomfort, irritation, pressure, or redness in this area, contact your cochlear implant professional.

- Although electromagnetic emissions of the Slim HPs have been confirmed within safe limits, some other devices could be sensitive to such emissions. If you notice other devices do not behave as expected when your Slim HPs are nearby, separate the devices.*
- The magnetic field created by the headpiece magnet may alter the programming of magnetic shunts. Users of magnetic shunts should consider off-the-shelf retention methods.*
- In the event that you experience any issues with your product, please contact your cochlear implant professional or the manufacturer. Do not attempt to service or modify the Slim HPs. Doing so may compromise system performance and will void the manufacturer's warranty. Products should be serviced only at Advanced Bionics and damaged products should be returned to Advanced Bionics.*

Undesirable Side Effects

Undesirable side effects of your Slim HPs may include skin irritation and discomfort from pressure on the ear or head, device overheating, or overly loud sounds.

Product Description

The Slim HP is a headpiece with an integrated cable designed for use with a sound processor worn on the ear. The Slim HP is available with cable lengths of 3.5 inches/9cm, 4.25 inches/11cm, and 5.5 inches/14cm.


*Sound Processor
Connector*

*Integrated
Cable*

Slim HP with Color Cap

The Slim HP Mic is a headpiece with an integrated 12 inch/30cm cable and a microphone designed for use with a sound processor worn off the ear (for example, in a clip).


*Slim HP Mic with
Color Cap*

Microphone

*Sound Processor
Connector*

Integrated Cable

The Slim HP AquaMic™ headpiece is a waterproof headpiece with an integrated cable and a microphone designed for use with waterproof solutions (e.g., Neptune, AquaCase, or M Waterproof Battery). The Slim HP AquaMic is available with cable lengths of 12 inches/30cm, 18 inches/46cm, 24 inches/61cm, and 42 inches/107cm.

The Slim HPs are compatible with the following accessory products: Slim HP Magnets, Slim HP Color Caps, and the Slim HP Color Cap Removal Tool.

The Slim HPs contain a magnet well designed to accommodate the appropriate Slim HP magnet with appropriate strength to ensure retention between the cochlear implant and the headpiece. The Standard Magnets are for use by recipients of CII, HiRes 90k, HiRes 90k Advantage, and HiRes Ultra implants. The 3D Magnets are for use by recipients of HiRes Ultra 3D implants.

The Slim HP Color Cap snaps in place on top of the headpiece securely covering the magnet well and the microphone. There are two sizes of Slim HP Color Caps: Standard and Large. The Standard Color Cap is for use with the smallest Slim HP Magnets: Standard Magnet sizes 1

through 5 and 3D Magnet sizes 1 through 4. The Large Color Cap is for use with the largest Slim HP Magnets: Standard Magnet sizes 6 and 7 and 3D Magnet sizes 5, 6, and 7.

INDICATIONS FOR USE

The Slim HP, Slim HP Mic, Slim HP AquaMic, and Slim HP Color Caps are external components of the HiResolution Bionic Ear System, which is intended to restore a level of auditory sensation to individuals with severe-to-profound sensorineural hearing loss via electrical stimulation of the auditory nerve.

Adults

- *18 years of age or older.*
- *Severe-to-profound bilateral sensorineural hearing loss or severe-to-profound unilateral hearing loss.*
- *Post-lingual onset of severe or profound hearing loss.*
- *Limited benefit from appropriately fitted hearing aids, defined as scoring 50% or less on a test of open-set sentence recognition (HINT Sentences).*

Children

- *12 months through 17 years of age.*
- *Severe-to-profound bilateral sensorineural deafness or severe-to-profound unilateral hearing loss.*
- *Use of appropriately fitted hearing aids for at least 6 months in children 2 through 17 years of age, or at least 3 months in children 12 through 23 months of age. The minimum duration of hearing aid use is waived if x-rays indicate ossification of the cochlea.*

- *Little or no benefit from appropriately fitted hearing aids. In younger children (< 4 years of age), lack of benefit is defined as a failure to reach developmentally appropriate auditory milestones (such as spontaneous response to name in quiet or to environmental sounds) measured using the Infant-Toddler Meaningful Auditory Integration Scale or Meaningful Auditory Integration Scale or \leq 20% correct on a simple open-set word recognition test (Multisyllabic Lexical Neighborhood Test) administered using monitored live voice (70 dB SPL). In older children (\geq 4 years of age), lack of hearing aid benefit is defined as scoring \leq 12% on a difficult open-set word recognition test (Phonetically Balanced-Kindergarten Test) or \leq 30% on an open-set sentence test (Hearing In Noise Test for Children) administered using recorded materials in the sound field (70 dB SPL).*

Intended Purpose

The Slim HP headpieces are part of the external components of an auditory active-implantable system, the HiResolution Bionic Ear system. The HiResolution Bionic Ear system is intended to provide auditory sensation via electrical stimulation of the auditory nerve for individuals with severe to profound bilateral or unilateral sensorineural hearing loss. Severe hearing loss is defined as audiometric thresholds greater than or equal to 70 dB HL, but less than 90 dB HL. Profound hearing loss is defined as audiometric thresholds greater than or equal to 90 dB HL. The external components work together with the implant to bypass the damaged part of the inner ear and convert sound picked up by the microphone or streamed via wireless communication into electrical signals that are used by the cochlear implant to enable hearing.

The Slim HP is a headpiece with an integrated cable designed for use with a sound processor on the ear and is intended to transfer sound signals and power between Advanced Bionics' cochlear implant and sound processor.

The Slim HP Mic headpieces are part of the external components of an auditory active-implantable system, the HiResolution Bionic Ear system. The HiResolution Bionic Ear system is intended to provide auditory sensation via electrical stimulation of the auditory nerve for individuals with severe to profound bilateral or unilateral sensorineural hearing loss. Severe hearing loss is defined as audiometric thresholds greater than or equal to 70 dB HL, but less than 90 dB HL. Profound hearing loss is defined as audiometric thresholds greater than or equal to 90 dB HL. The external components work together with the implant to bypass the damaged part of the inner ear and convert sound picked up by the microphone or streamed via wireless communication into electrical signals that are used by the cochlear implant to enable hearing.

The Slim HP Mic is a headpiece with an integrated cable and a microphone designed for use with a sound processor worn off the ear. It is intended to transfer sound signals and power between Advanced Bionics' cochlear implant and sound processor. The embedded microphone, if enabled, allows the headpiece to also function as a sound input source.

The Slim HP AquaMic™ headpieces are part of the external components of an auditory active-implantable system, the HiResolution Bionic Ear system. The HiResolution Bionic Ear system is intended to provide auditory sensation via electrical stimulation of the auditory nerve for individuals

with severe to profound bilateral or unilateral sensorineural hearing loss. Severe hearing loss is defined as audiometric thresholds greater than or equal to 70 dB HL, but less than 90 dB HL. Profound hearing loss is defined as audiometric thresholds greater than or equal to 90 dB HL. The external components work together with the implant to bypass the damaged part of the inner ear and convert sound picked up by the microphone or streamed via wireless communication into electrical signals that are used by the cochlear implant to enable hearing.

The Slim HP AquaMic™ is a headpiece with an integrated cable and a microphone designed for use in water environments. It is intended to transfer sound signals and power between Advanced Bionics' cochlear implant and sound processor. The embedded microphone, if enabled, allows the headpiece to also function as a sound input source.

Intended Users

The intended users of the Slim HPs are recipients of the Neptune, Naída™ CI Q, Naída™ CI M, or Sky CI™ M sound processors, their caregivers if applicable, and hearing care professionals.

The recipient, or their caregiver, should, at a minimum, be capable of changing the battery, connecting the earhook, placing and removing the sound processor from their ear, and placing and removing the headpiece from the implant site.

The hearing care professional should be trained in the use and fitting of the cochlear implant system.

Intended Use Environment

The intended use environments of the Slim HPs are healthcare and daily living environments. The Slim HP, Slim HP Mic, Slim HP AquaMic, and Slim HP Color Caps are designed to be worn daily on the skin over the implant during the user's normal waking hours.

Compatibility

The Slim HPs are compatible with Neptune, Naída CI Q, Naída CI M, or Sky CI M sound processors from Advanced Bionics.

Operating Instructions

Assembly Instructions

To operate your headpiece with your cochlear implant system you will need to install a Slim HP Magnet and attach a Slim HP Color Cap, as follows:

- 1.** Insert the Slim HP Magnet into the magnet well of the headpiece.
 - a.** For recipients using the Standard Magnet, no specific magnet orientation is required.
 - b.** For recipients using the 3D Magnet, the magnet should be inserted so that the arrow on the magnet points up when the headpiece is on the head.
- 2.** Attach a color cap (See “Color Caps” section).
- 3.** Attach the headpiece and cable to your sound processor (See “Cables” section).

4. Then place the bottom side of the headpiece over the internal magnet of your cochlear implant.

NOTE: For recipients of HiRes Ultra 3D implants, if the headpiece repels or pushes away when brought over the internal device, rotate the headpiece back and forth as you move the headpiece closer towards the implant. Once attraction is felt, adjust the headpiece to the original orientation and place on head.

NOTE: For recipients of HiRes Ultra 3D implants, when in a magnetic field or when the headpiece is removed, you may feel the internal components of the magnet vibrate and/or rotate. This movement is expected as the magnets are designed to rotate in the presence of strong magnetic fields. This movement may be noticeable and is not expected to be uncomfortable.

Cables

Slim HP, Slim HP Mic, and Slim HP AquaMic have cables that are integrated into the headpiece.

Attaching the Slim HP cable to the sound processor

To attach the Slim HP cable to the sound processor, hold the headpiece cable by the cable's strain relief. Align the connector on the cable with the connector on the sound processor and gently push the cable into the sound processor until it snaps into place.

Removing the Slim HP cable from the sound processor

Remove the battery from the processor before removing the cable from the processor. To remove the cable from the sound processor, always hold the cable's strain relief and gently pull it away.

NOTE: To maximize cable life, the cable should only be unplugged from the sound processor when changing wearing configuration or when replacing the cable.

Color Caps

Attaching the Slim HP Color Cap

Align the Color Cap over the Slim HP and snap it into place.

Removing the Slim HP Color Cap

The arrow on the underside of the Slim HPs indicates where the Slim HP Color Cap removal tool can be inserted to remove the Color Cap from the Slim HP. Insert the tool as indicated by the arrow to pop the Color Cap off of the Slim HP.

Magnets

It is important for your headpiece to have the correct magnet strength and orientation to avoid discomfort or retention issues. Contact your cochlear implant professional if there are any concerns regarding magnet retention. If deemed appropriate, your cochlear implant professional may replace the headpiece magnet with a stronger or weaker magnet or adjust magnet orientation. If needed, you may use off-the-shelf solutions to improve retention.

For use of Slim HP AquaMic in water

If a swim cap, headband, or other device used in water environments is placed over the headpiece microphone, some degradation in sound quality may be perceived. If the headpiece is worn following submersion in water, sound quality may be affected. Dry the headpiece fully to restore typical system sound quality. It is recommended that you replace your headpiece after water use with a dry headpiece to ensure restoration of typical system sound quality.

CARE AND MAINTENANCE

Storing the Slim HPs

When not in use, store Slim HPs in the equipment case provided.

Cleaning and Maintenance

Periodically, clean your headpiece using a soft cloth. Do not use solvents. Replace the color cap every 90 days, or sooner, if a change in sound quality is noted. Following use of the Slim HP AquaMic in Water:

- Rinse Slim HP AquaMic with clean water.

- Remove the color cap.
- Wipe or blot the headpiece upper surface with a soft cloth.
- Place in an appropriate drying system to dry completely.

Unaided listeners may periodically assess headpiece sound quality using the listening check function. For instructions, consult your appropriate sound processor user guide.

Recommended Operating and Storage Temperature Ranges

CONDITION	MINIMUM	MAXIMUM
Operating Temperature	0°C (32°F)	45°C (115°F)
Storage Temperature	-20°C (-13°F)	55°C (131°F)
Relative Humidity	0%	95%

Expected Life

The expected life of the Slim HP, Slim HP Mic, and Slim HP AquaMic is 5 years. The expected life of the Slim HP Color Caps is 90 days.

Ingress Protection (IP) Ratings

The following product configuration has an IP rating of 57 (protection against dust; protection against failure due to one-time immersion for 30 minutes to a depth up to 1 meter and after drying overnight in a drying system):

- **The Slim HP or Slim HP Mic and a Naída CI Q or Naída CI M sound processor.**

The following product configuration has an IP rating of 68 (complete protection against dust penetration; protection against failure due to continuous immersion in water up to 3 meters):

- **The Slim HP AquaMic with a Neptune sound processor, AquaCase, or M Waterproof Battery.**

DISPOSAL AND HAZARDOUS MATERIALS

Dispose of the Slim HPs in accordance with applicable national and local regulations. Hazardous materials: not applicable to this product.

AVAILABLE PRODUCT OPTIONS AND MODEL NUMBERS

Slim HP

PRODUCT OPTION	MODEL NUMBER
Beige 3.5"/9cm Cable	CI-5320-001
Beige 4.25"/11cm Cable	CI-5320-006
Beige 5.5"/14cm Cable	CI-5320-011
Black 3.5"/9cm Cable	CI-5320-002
Black 4.25"/11cm Cable	CI-5320-007
Black 5.5"/14cm Cable	CI-5320-012
White 3.5"/9cm Cable	CI-5320-005

PRODUCT OPTION	MODEL NUMBER
White 4.25"/11cm Cable	CI-5320-010
Gray 3.5" /9cm Cable	CI-5320-004
Gray 4.25"/11cm Cable	CI-5320-009
Brown 3.5"/9cm Cable	CI-5320-003
Brown 4.25"/11cm Cable	CI-5320-008

Slim HP Mic

PRODUCT OPTION	MODEL NUMBER
Black 12"/30cm Cable	CI-5322-001

Slim HP AquaMic

PRODUCT OPTION	MODEL NUMBER
Beige 12"/30cm Cable	CI-5321-001
Beige 18"/46cm Cable	CI-5321-005
Beige 24"/61cm Cable	CI-5321-009
Beige 42"/107cm Cable	CI-5321-013
Black 12"/30cm Cable	CI-5321-002

PRODUCT OPTION	MODEL NUMBER
Black 18"/46cm Cable	CI-5321-006
Black 24"/61cm Cable	CI-5321-010
Black 42"/107cm Cable	CI-5321-014
Gray 12"/30cm Cable	CI-5321-004
Gray 18"/46cm Cable	CI-5321-008
Gray 24"/61cm Cable	CI-5321-012
Gray 42"/107cm Cable	CI-5321-016
Brown 12"/30cm Cable	CI-5321-003
Brown 18"/46cm Cable	CI-5321-007
Brown 24"/61cm Cable	CI-5321-011
Brown 42"/107cm Cable	CI-5321-015

Slim HP Color Caps

PRODUCT OPTION	MODEL NUMBER
Standard Sand Beige	CI-7131-004
Standard Chestnut	CI-7131-005
Standard Silver Grey	CI-7131-006

PRODUCT OPTION	MODEL NUMBER
Standard Velvet Black	CI-7131-007
Standard Alpine White	CI-7131-010
Standard Lava Red	CI-7131-001
Standard Blue Ocean	CI-7131-002
Standard Majesty Purple	CI-7131-003
Standard Caribbean Pirate	CI-7131-008
Standard Precious Pink	CI-7131-009
Standard Resin Black	CI-7131-011
Large Sand Beige	CI-7132-001
Large Chestnut	CI-7132-002
Large Silver Grey	CI-7132-003
Large Velvet Black	CI-7132-004
Large Resin Black	CI-7132-005
Color Cap Removal Tool	CI-7526

CLINICAL BENEFITS

The Intended Clinical Benefit of the Slim HP, Slim HP Mic, Slim HP AquaMic, and Slim HP Color Caps as parts of the HiResolution Bionic Ear System is to:

- *provide useful hearing to individuals with severe-to-profound hearing loss via electrical stimulation of the auditory nerve.*
- *provide a combined type of stimulation, electric plus acoustic, when sound processors are enabled with the Acoustic Earhook.*

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

***A summary of safety and clinical performance (SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI:- Slim HP - 08400944CI5320XQ, Slim HP AquaMic - 08400944CI5321XS, Slim HP Mic - 08400944CI5322XU
<https://ec.europa.eu/tools/eudamed>***

Users of the Slim HPs in the European Union should report any serious incident to their local competent authority and to Advanced Bionics.

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©2021 Advanced Bionics LLC. All Rights Reserved.
Advanced Bionics LLC