INSTRUCTIONS FOR USE  
**HiFocus™ Mid-Scala Reusable Surgical Tools**  
For more information contact:

 Advanced Bionics LLC  
28515 Westinghouse Place  
Valencia, CA 91355, USA  
Tel: + 1 661-362-1400

**EC|REP** Advanced Bionics GmbH  
Feodor-Lynen-Strasse 35  
D-30625 Hannover  
Tel: +49 511 5248 7522

[www.advancedbionics.com](http://www.advancedbionics.com)

## Contents

Symbols Table	2
Purpose and Scope of IFU	2
Device Description	3
Intended Purpose	3
Target/Patient Population	3
Intended Users	3
Information for Use and Required Training	3
Intended Use Environment	4
Indications for Use	4
Limitations and Contraindications	4
Precautions, Cautions, and Warnings	5
Possible Adverse Events	6
MRI Safety Information	6
Disposal and Hazardous Materials	6
Expected Life	7
Clinical Benefits	7
Summary of Safety and Clinical Performance	8
Instructions for Use	8
Care, Cleaning, and Maintenance	9
Available Product Options and Model Numbers	14

## Symbols Table

SYMBOL	EXPLANATORY TEXT
	European Community Mark of Conformity. Authorized to affix the CE Mark in 2022.
	Date of Manufacture
	Manufacturer
	Model Number
	Lot Number
	European Authorized Representative
	Unique Device Identifier
	Medical Device
	Use-by Date
	Suitable temperature range for transport and storage
	See Instructions for Use
	Caution (Refer to “Cautions and Warnings” section for further information)
	MR Unsafe

## Purpose and Scope of IFU

This IFU is designed to help cochlear implant hearing care professionals trained in the implantation of cochlear implants to understand the use and care of the HiFocus™ Mid-Scala Reusable Surgical Tools Kit.

## **Device Description**

The HiFocus™ Mid-Scala Electrode Instrument Kit (CI-4508) is a surgical kit containing two sets of reusable tools and gauges for insertion of the HiFocus Mid-Scala Electrode: Two (2) identical cochleostomy gauges (CI-4347) and two (2) identical Electrode Claw Tools (CI-4254).

The HiFocus™ Mid-Scala Cochleostomy Gauge (CI-4347) is a non-sterile reusable gauge labeled “0.8 mm” that can be used to verify the minimum recommended cochleostomy or round window size. The tools must be sterilized prior to use.

The HiFocus™ Mid-Scala Claw Tool (CI-4254) is a non-sterile tool that can be used to stabilize the electrode during and/ or following insertion or to push the array in during freehand insertion. The tools must be sterilized prior to use.

## **Intended Purpose**

The HiFocus™ Mid-Scala Electrode Instrument Kit (CI-4508) are intended to provide the tools and gauges needed to assist with insertion of the HiFocus Mid-Scala Electrode.

The HiFocus™ Mid-Scala Cochleostomy Gauge (CI-4347) is intended to verify the minimum recommended cochleostomy size.

The HiFocus™ Mid-Scala Claw Tool (CI-4254) is intended to stabilize the electrode during and/or following insertion.

Both tools are tissue contacting and are made from 304 stainless steel.

## **Target/Patient Population**

The target population for the HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool are cochlear implant candidates and recipients.

## **Intended Users**

The intended users of the HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool are ENT surgeons trained in the implantation of cochlear implants.

## **Information for Use and Required Training**

A Surgeon’s Manual and a video describing the surgical procedure and insertion of the electrode are provided to all ENT surgeons prior to implantation. ENT surgeons must be well versed in mastoid surgery and the facial recess approach to the round window. Advanced Bionics conducts periodic training courses on the recommended surgical procedure to implant HiRes™ Ultra and HiRes™ Ultra 3D. ENT surgeons who implant adults and/or pediatrics must receive training before implantation. ENT surgeons should work with an audiology professional who has been trained fully on

the proper fitting and adjustment of the system. Failure to obtain the appropriate training might result in higher incidence of surgical and medical complications.

## **Intended Use Environment**

The intended use environments of the HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool are sterile operating rooms.

## **Indications for Use**

The HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool are surgical tool are used to assist in insertion of the Mid-Scala Electrode process during implantation of the HiResolution Bionic Ear System.

## **Limitations and Contraindications**

As the HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool is a surgical tool used for the implantation of the HiResolution Bionic Ear System, the following contraindications stated for the HiResolution Bionic Ear System are applicable: Deafness due to lesions of the acoustic nerve or central auditory pathway; cochlear ossification that prevents electrode insertion; absence of cochlear development; tympanic membrane perforations associated with recurrent middle ear infections.

## Precautions, Cautions, and Warnings

Bacterial *meningitis* has been reported in users of the system and other cochlear implants, especially in children under the age of 5. The cause of meningitis in these cases has not been established. A small percentage of deaf patients may have congenital abnormalities of the cochlea (inner ear) which predispose them to meningitis even prior to implantation. Patients who become deaf as a result of meningitis are also at increased risk of subsequent episodes of meningitis compared to the general population. Other predisposing factors may include young age (<5 years), otitis media, immunodeficiency, or surgical technique. The cochlear implant, because it is a foreign body, may act as a nidus for infection when patients have bacterial illnesses.

The incidence rate, although low, appears to be higher than the age- adjusted rate for the general population. The fatality rate as a result of meningitis also appears to be higher. Adequate epidemiological data are not available to determine whether the incidence and fatality rates are, in fact, definitively different from the general population, whether there are special risk factors in the cochlear implant population, or whether different cochlear implant models pose different risks.

Adults and parents of children who are considering a cochlear implant or who have received cochlear implants should be advised of the risk of meningitis. They should also be informed of the availability of vaccines that have been shown to substantially reduce the incidence of meningitis in the general population resulting from the organisms that commonly cause bacterial meningitis (*Streptococcus pneumoniae*, *Haemophilus influenzae*, *Meningococcus*). National health agencies frequently provide updated information on the safety and utility of specific vaccines and offer recommendations reflecting local or regional conditions. Physicians or patients should refer to the applicable authorities for this information. These vaccines can be administered by pediatricians, primary care/family physicians, and infectious disease specialists.

Adults and parents of children who have received cochlear implants should be counseled on the symptoms of meningitis, the need to seek immediate medical care if any symptoms appear, and the need to advise the treating physicians of the presence of the cochlear implant and of the possibility of increased risk of meningitis associated with implant. They should also be counseled to obtain medical care at the first signs of otitis media.

**Electrosurgery:** Electrosurgical instruments and RF Ablation instruments are capable of producing radio-frequency voltages of such magnitude that a direct coupling might occur between the cautery tip and the electrode. Induced currents could cause damage to the cochlear tissues or permanent damage to the implant. Prior to placing the implant in the patient, ensure that electrosurgery has been discontinued.

## Possible Adverse Events

The following risks associated with cochlear implantation and ear surgery also can occur:

- Implant patients incur the normal risks of surgery and general anesthesia.
- Major ear surgery may result in numbness, swelling or discomfort about the ear, disturbance of taste or balance, or neck pain. If these events occur, they are usually temporary and subside within a few weeks of surgery.
- Rarely, cochlear implantation may cause a leak of the inner ear fluid, which may result in meningitis.
- During the surgery, it is a rare possibility that the facial nerve could be injured resulting in a temporary or permanent weakening or full paralysis on the same side of the face as the implant.
- During the surgery, there is a rare possibility that cerebrospinal fluid leakage or perilymph fluid leakage could occur.
- As a result of the surgery, it is possible that dizziness, tinnitus, or vertigo may result. If these events occur, they are usually temporary and subside over time.
- There is a possibility that the electrode or device may migrate due to improper insertion, requiring additional medical treatment or removal of the internal device to address any resulting injury.

*Users of the HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool should report any serious incident to their competent authority and to Advanced Bionics.*

## MR Safety Information

The components of the HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool are MR Unsafe. Insertion of the HiFocus™ Mid-Scala Electrode is contraindicated in an MR environment.

## Disposal and Hazardous Materials

Dispose of the HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool and accessories in accordance with applicable national and local regulations for disposing of biohazardous materials.

## **Expected Life**

The expected life of the HiFocus™ Mid-Scala Electrode Instrument Kit is no more than 20 cycles of reprocessing. If using a pre-vacuum sterilization cycle of 134°C (273°F) at 18 minutes, the expected life is no more than 3 cycles of reprocessing.

The expected life of the HiFocus™ Mid-Scala Cochleostomy Gauge is no more than 20 cycles of reprocessing. If using a pre-vacuum sterilization cycle of 134°C (273°F) at 18 minutes, the expected life is no more than 3 cycles of reprocessing.

The expected life of the HiFocus™ Mid-Scala Claw Tool is no more than 20 cycles of reprocessing. If using a pre-vacuum sterilization cycle of 134°C (273°F) at 18 minutes, the expected life is no more than 3 cycles of reprocessing.

## **Clinical Benefits**

The Intended Clinical Benefit of the HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool as part of the HiResolution Bionic Ear System is to assist with the implantation of the Mid-Scala Electrode.

## Summary of Safety and Clinical Performance

A summary of safety and clinical performance (SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI:

### CI-4508 - HiFocus™ Mid-Scala Electrode Instrument Kit

- 08400944CI4508Y5

### CI-4347 - HiFocus™ Mid-Scala Cochleostomy Gauge

- 08400944CI42347TF

### CI-4254 - HiFocus™ Mid-Scala Claw Tool

- 08400944CI4254XV

<https://ec.europa.eu/tools/eudamed>

Users of the HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool in the European Union should report any serious incident to their local competent authority and to Advanced Bionics.

## Instructions for Use

The HiFocus™ Mid-Scala Electrode Instrument Kit (CI-4508) includes the following tools and gauges for insertion of the HiFocus™ Mid-Scala Electrode:

- Two (2) Cochleostomy Gauges (CI-4347) labeled “0.8 mm”
- Two (2) Electrode Claw Tools (CI-4254)

The cochleostomy gauge can be used to verify the minimum recommended cochleostomy or round window size.

The Claw Tool can be used to stabilize the electrode during and/or following insertion or to push the array in during freehand insertion. The Claw Tool may help prevent buckling/movement of the electrode.

# Care, Cleaning, and Maintenance

## Warnings and Cautions

- All tools are supplied unsterilized, therefore, all tools must be cleaned, disinfected and sterilized prior to use.
- Sterilization is always required regardless of which of the disinfection method steps is used.
- Personal Protective Equipment should be worn when handling or working with contaminated or potentially contaminated materials, devices, equipment and cleaning agents.
- Process temperatures should not exceed 145°C (293°F) to prevent damage to the tools.
- Do not use harsh chemicals such as chlorine, caustic soda, and organic or ammoniated acids or solvents (such as acetone), as these are likely to damage the tools.
- The cleaning/disinfection agent manufacturer's recommendation for the respective specifically applied agent needs to be followed.
- Thermal disinfection may only be performed in connection with automated cleaning, applying the disinfection step after rinsing and prior to drying steps.

## Reprocessing Guidelines

The reprocessing of the tools is split into five stages: initial treatment at the point of use, preparation before cleaning, cleaning, disinfection, sterilization. Prior to each use thereafter, be sure to follow the five stages outlined in these reprocessing guidelines.

### I: Initial Treatment at the Point of Use

- At the point of use and after use, submerge tools in water; this will prevent any contaminate drying on the tools.
- It is recommended that tools be reprocessed as soon as possible after use.

### II: Preparation Before Cleaning

- Rinse the tools under warm running tap water<sup>1</sup> until visible soil is removed from the tools.

- Equipment for Cleaning Methods: soft bristle brush<sup>2</sup>, ultrasonic unit<sup>3,4</sup>, lint-free cloth, alkaline detergent with pH≥10 -- nonabrasive, low-foaming<sup>5</sup>.

### III: Cleaning

- For cleaning, either the manual or automatic method can be used.

#### Manual Processing Steps:

1. Prepare alkaline detergent per manufacturer's recommendation using warm tap water<sup>1</sup>.
2. Fully submerge the tools in alkaline detergent solution and soak for three minutes.
3. After three minutes of soaking, brush the tools with the soft bristle brush for one minute to remove any debris or soil.
4. Remove tools from the detergent solution and rinse under warm running tap water<sup>1</sup> for three minutes.
5. Prepare alkaline detergent per the manufacturer's recommendation using warm tap water<sup>1</sup> in ultrasonic unit.
6. Fully submerge the tools in the ultrasonic unit and sonicate for 15 minutes.
7. Important: If reprocessing multiple tools, ensure the tools are not in contact with each other during sonication.
8. Remove the tools from the ultrasonic unit at the completion of sonication.
9. Rinse the tools with Reverse Osmosis (RO) or Deionization (DI) water for one minute.
10. Dry the tools with lint-free cloth and then visually inspect tools for visual soil. If soil is present, repeat the cleaning process.

#### Automatic Processing Steps:

1. Rinse the tools under warm running tap water<sup>1</sup> until no visible soil is present.
2. A soft bristle brush or syringe may be used to aid in the rinsing.
3. Place the tool in the washer-disinfector according to the following parameters:

### Minimum Washer-Disinfector Cycle Parameters for Non-EU Countries

Phase	Recirculation Time (Minutes)	Water	Temperature	Detergent
Pre-Wash	2:00	Cold Tap Water	65.5°C (150°F)	None
Wash	2:00		65.5°C (150°F)	Alkaline Detergent
Rinse	1:00		65.5°C (150°F)	None
For Thermal Disinfection, please see Step IV				
Drying	15:00		90°C (194°F)	None

4. Dry the tools with a lint-free cloth and then visually inspect the tools for visual soil. If soil is present, repeat the cleaning process.

#### IV: Disinfection

- For disinfection, either the chemical or thermal method can be used.

##### Chemical Processing Steps:

1. Equipment: 10 mL syringes, lint-free cloth, Cidex® 2.4% Glutaraldehyde<sup>6</sup> and Cidex 2 test strips<sup>7</sup>.
2. Prepare the disinfectant per the manufacturer's recommendation using test strips to determine concentration.
3. Fully submerge the tools in the disinfectant for 60 minutes at 25°C ± 2°C (77°F ± 3.6°F).
4. Remove any air bubbles from tools surface using a lint-free cloth or a syringe.
5. After 60 minutes, remove the tools from the disinfectant solution and place in a DI water bath.
6. Submerge the tools in the DI water bath for 30 minutes, use a 10 mL (0.3 oz.) syringe to force 10 mL (0.3 oz.) of DI water into cracks and grooves for a total of six times.
7. Repeat processing step 5 one additional time with a fresh bath of DI water.

## Thermal Processing Steps:

1. Equipment: Washer-Disinfector<sup>8,9</sup>
2. Thermal disinfection may only be performed in connection with automated cleaning, applying the disinfection step after rinsing and prior to drying steps.
3. For automated cleaning, use the following parameters:

Phase	Minimum Water Temperature	Minimum Recirculation Time
Thermal	90°C (194°F)	05:00 minutes

4. IMPORTANT: After cleaning and disinfection, use appropriate packaging compliant to ISO 11607-1: 2019, as applicable, before sterilization.

## V: Sterilization

Equipment: Prevacuum sterilizer<sup>10</sup>, Gravity Displacement Sterilizer<sup>11</sup>, CSR wrap<sup>12</sup>, Self-seal Pouches<sup>13,14</sup>

### Processing Steps:

Wrap tools in a double barrier (pouches or sterilization wrap using the sequential envelope folding technique), and use the following parameters (in accordance with local guidelines):

<b>Sterilizer Cycle Parameters for EU Countries</b>			
Sterilizer Type	Temperature	Dwell Time	Drying
Pre-Vacuum	134°C (273°F)	3 minutes	30 minutes
	134°C (273°F)	18 minutes	30 minutes

<b>Sterilizer Cycle Parameters for Non-EU Countries</b>			
Sterilizer Type	Temperature	Dwell Time	Drying
Pre-Vacuum	132°C (270°F)	4 minutes	30 minutes
	134°C (273°F)	3 minutes	30 minutes
	134°C (273°F)	18 minutes	30 minutes
Gravity Displacement	134°C (273°F)	15 minutes	30 minutes

**IMPORTANT:** When sterilizing multiple tools, ensure that the sterilizer's maximum load is not exceeded.

## INSPECTION AND MAINTENANCE

- Carefully examine each tool to ensure that all visible contamination has been removed. If contamination is noted repeat the cleaning/disinfection process.
- Visually inspect the tool(s) and if any of the following is noted discard the tool(s) and contact Advanced Bionics as the function or use of the tool(s) may be compromised:
  - Corrosion
  - Discoloration
  - Pitting or cracks
  - Markings are not legible
- Examine the tool(s) to ensure its functional. If found as not functional, discard the tool(s) and contact Advanced Bionics.

## STORAGE RECOMMENDATIONS

Before sterilization, use appropriate packaging compliant to ISO 11607-1: 2019, as applicable.

It is recommended that you follow the storage process established by the processing facility.

---

The instructions provided above have been validated by Advanced Bionics as being capable of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the processing as actually performed using standards, equipment, materials and personnel in the processing facility achieves the desired result. This requires validation and routine monitoring of the process.

Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Provisions of specific countries may deviate from the presented process and shall additionally be taken into account during reprocessing. This guide is provided in fulfillment of EN ISO 17664:2021, EN ISO 15883 -1: 2009, EN ISO 11607-1: 2019 and ANSI/AAMI ST79: 2017. The processor's large steam sterilizers are to have feed water compliant with EN 285:2015. The compliance of the cited standards by the processor is a prerequisite to the recommended process.

The following materials were used to validate the referenced processes:

1. Water Temperature Ranges used in the Validation: Warm water (27.1°C – 31.4°C) and Cold Water (21.0°C – 27.0°C).
2. Model M-16 Soft Bristle Brush.
3. Model Branson® CPXH Ultrasonic Bath.
4. Model Branson® M8800 Ultrasonic Bath.
5. Neodisher Mediclean Forte Detergent (concentration 1/4 oz./g (2mL/L).
6. Cidex® 2.4% Glutaraldehyde.
7. Cidex® 2.4 Test Strips.
8. Model SD880 Reliance Vision™ Single-Chamber Washer/Disinfector.
9. Model Belimed WD290 Washer-Disinfector.
10. Model Steris Amsco Lab 250 Prevacuum Sterilizer.
11. Model Steris Amsco Lab 250 Gravity Displacement Sterilizer.
12. Model 4040 Cardinal Health Convertors BioShield Regular, CSR Wrap.
13. Model 92308 Cardinal Health Self-seal Pouch.
14. Model 92510 Cardinal Health Self-seal Pouch.

### **Available Product Options and Model Numbers**

The HiFocus™ Mid-Scala Electrode Instrument Kit, HiFocus™ Mid-Scala Cochleostomy Gauge, and HiFocus™ Mid-Scala Claw Tool is compatible with the following implant types from Advanced Bionics:

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