INSTRUCTIONS FOR USE
HiResolution™ Bionic Ear System, HiRes™ Ultra 3D HiFocus™ MS & SlimJ Electrodes

For more information contact:

 Advanced Bionics LLC
 28515 Westinghouse Place
 Valencia, CA 91355 USA
 +1 661-362-1400

 Advanced Bionics GmbH
 Feodor-Lynen-Strasse 35
 D-30625 Hannover
 +49 511 5248 7522

 0123
 2021
 ifu.advancedbionics.com

Contents

SYMBOLS TABLE 1

PURPOSE AND SCOPE OF IFU 1

DEVICE DESCRIPTION 1

INTENDED PURPOSE 2

INTENDED USERS 2

INTENDED USE ENVIRONMENT 2

TARGET/PATIENT POPULATION 2

INDICATIONS FOR USE 2

CONTRAINDICATIONS 2

PRECAUTIONS, CAUTIONS, AND WARNINGS 2

MRI SAFETY INFORMATION 3

PRECAUTIONS 3

EXPECTED LIFE 4

CLINICAL BENEFITS 4

CLINICAL STUDIES 4

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE 4

POSSIBLE ADVERSE EVENTS 4

PATIENT COUNSELING INFORMATION 4

INFORMATION FOR USE AND REQUIRED TRAINING 5

SPECIAL HANDLING 5

SYMBOLS TABLE

| SYMBOL | EXPLANATORY TEXT |
|---|--|
|  | European Community Mark of Conformity. Authorized to affix the CE Mark in 2021 |
|  | Date of Manufacture |
|  | Manufacturer |
|  | Model number |
|  | Serial number |
|  | European authorized representative |
|  | Unique Device Identifier |
|  | Medical Device |
|  | Use-by Date |
|  | Sterile Indication and Method |
|  | Double barrier sterile pack sterilized with ethylene oxide |
|  | Sterile Lot Barcode |
|  | Single Use / Do Not Re-Use |
|  | Suitable temperature range for transport and storage |
|  | Suitable pressure range for transport and storage |
|  | See Instructions for Use |
|  | Fragile |
|  | Caution (refer to "Cautions & Warnings" section for further information) |
|  | MR Conditional |

PURPOSE AND SCOPE OF IFU

This IFU is designed to help cochlear implant recipients, and their caregivers, if applicable, to understand the use of the HiRes™ Ultra 3D HiFocus™ Mid-Scala (MS) Electrode and HiRes™ Ultra 3D HiFocus™ SlimJ Electrode cochlear implants (HiRes™ Ultra 3D). It is meant to be used by Ear, Nose and Throat (ENT) surgeons and hearing care professionals (audiologists) to provide system level instructions.

DEVICE DESCRIPTION

The HiResolution™ Bionic Ear System consists of internal and external components. The internal components include the HiRes™ Ultra 3D receiver and either the HiFocus™ SlimJ electrode or the HiFocus™ MS electrode array that are implanted surgically under the skin behind the ear. Both the HiFocus™ SlimJ and the HiFocus™ MS electrode array have 16 contacts and are connected to the implant through the electrode lead. The main difference between the HiFocus™ MS electrode array and the HiFocus™ SlimJ electrode is the curved versus straight design. The receiver stimulator encapsulates the electronics into a titanium casing, the antenna coil allows for forward and backward telemetry between the implant and the external parts. The antenna coil also includes a magnet in order to retain the external headpiece.

The external components include a sound processor (body-worn or ear-level), a headpiece, and a cable. The system converts sound into electrical energy that activates the auditory nerve. The auditory nerve then sends information to the brain, where it is interpreted as sound.

INTENDED PURPOSE

The HiRes™ Ultra 3D HiFocus™ MS Electrode cochlear implant is an auditory active implantable device in the HiResolution™ Bionic Ear System. The HiResolution™ Bionic Ear System is intended to provide auditory sensation via electrical stimulation of the auditory nerve for individuals with severe to profound bilateral or unilateral sensorineural hearing loss. Severe hearing loss is defined as audiometric thresholds greater than or equal to 70dB HL, but less than 90 dB HL. Profound hearing loss is defined as audiometric thresholds greater than or equal to 90 dB HL. The HiRes™ Ultra 3D Ci HiFocus™ MS Electrode cochlear implant receives power and sound data over an inductively coupled link from the external sound processor system and converts the sound data into electrical stimulation which is delivered to the auditory nerve via the pre-curved electrode array to enable hearing. A self-aligning internal magnet allows the cochlear implant to be scanned at 1.5T and 3.0T at any orientation within the MRI scanner without bandaging and without magnet removal.

The HiRes™ Ultra 3D HiFocus™ SlimJ Electrode cochlear implant is an auditory active implantable device in the HiResolution™ Bionic Ear System. The HiResolution™ Bionic Ear System is intended to provide auditory sensation via electrical stimulation of the auditory nerve for individuals with severe to profound bilateral or unilateral sensorineural hearing loss. Severe hearing loss is defined as audiometric thresholds greater than or equal to 70dB HL, but less than 90 dB HL. Profound hearing loss is defined as audiometric thresholds greater than or equal to 90 dB HL. The HiRes™ Ultra 3D Ci HiFocus™ SlimJ Electrode cochlear implant receives power and sound data over an inductively coupled link from the external sound processor system and converts the sound data into electrical stimulation which is delivered to the auditory nerve via the lateral wall electrode array to enable hearing. A self-aligning internal magnet allows the cochlear implant to be scanned at 1.5T and 3.0T at any orientation within the MRI scanner without bandaging and without magnet removal.

INTENDED USERS

The intended users of the HiRes™ Ultra cochlear implants are recipients of Advanced Bionics cochlear implants, their caregivers, if applicable, ENT surgeons and hearing care professionals (audiologists).

The ENT surgeons and hearing care professionals (audiologists) should be trained in the implantation, use and fitting of the cochlear implant system. Please see section *Information For Use and Required Training* for training requirements.

INTENDED USE ENVIRONMENT

The intended use environments of the HiRes™ Ultra 3D cochlear implants are healthcare and daily living environments. The HiRes™ Ultra 3D cochlear implants are implanted under the skin behind the ear of the recipient.

TARGET/PATIENT POPULATION

Adults – 18 years of age or older
Children – 12 months through 17 years of age

INDICATIONS FOR USE

The HiRes™ Ultra 3D cochlear implant as an internal component of the HiResolution™ Bionic Ear System is intended to restore a level of auditory sensation to individuals with severe-to-profound sensorineural hearing loss via electrical stimulation of the auditory nerve.

Adults

- Severe-to-profound bilateral sensorineural hearing loss or severe-to-profound unilateral hearing loss.
- Postlingual onset of severe or profound hearing loss.
- Limited benefit from appropriately fitted hearing aids, defined as scoring 50% or less on a test of open-set sentence recognition (HINT Sentences).

Children

- Severe-to-profound bilateral sensorineural deafness or severe-to-profound unilateral hearing loss.
- Use of appropriately fitted hearing aids for at least 6 months in children 2 through 17 years of age, or at least 3 months in children 12 through 23 months of age. The minimum duration of hearing aid use is waived if x-rays indicate ossification of the cochlea.

- Little or no benefit from appropriately fitted hearing aids. In younger children (< 4 years of age), lack of benefit is defined as a failure to reach developmentally appropriate auditory milestones (such as spontaneous response to name in quiet or to environmental sounds) measured using the Infant-Toddler Meaningful Auditory Integration Scale or Meaningful Auditory Integration Scale or $\leq 20\%$ correct on a simple open-set word recognition test (Multisyllabic Lexical Neighborhood Test) administered using monitored live voice (70 dB SPL). In older children (≥ 4 years of age), lack of hearing aid benefit is defined as scoring $\leq 12\%$ on a difficult open-set word recognition test (Phonetically Balanced-Kindergarten Test) or $\leq 30\%$ on an open-set sentence test (Hearing In Noise Test for Children) administered using recorded materials in the soundfield (70 dB SPL).

CONTRAINDICATIONS

Deafness due to lesions of the acoustic nerve or central auditory pathway; active external or middle ear infections; cochlear ossification that prevents electrode insertion; absence of cochlear development; tympanic membrane perforations associated with recurrent middle ear infections.

PRECAUTIONS, CAUTIONS, AND WARNINGS

- Bacterial **meningitis** has been reported in users of the system and other cochlear implants, especially in children under the age of 5. The cause of meningitis in these cases has not been established. A small percentage of deaf patients may have congenital abnormalities of the cochlea (inner ear) which predispose them to meningitis even prior to implantation. Patients who become deaf as a result of meningitis are also at increased risk of subsequent episodes of meningitis compared to the general population. Other predisposing factors may include young age (<5 years), otitis media, immunodeficiency, or surgical technique. The cochlear implant, because it is a foreign body, may act as a nidus for infection when patients have bacterial illnesses.

The incidence rate, although low, appears to be higher than the age-adjusted rate for the general population. The fatality rate as a result of meningitis also appears to be higher. Adequate epidemiological data are not available to determine whether the incidence and fatality rates are, in fact, definitively different from the general population, whether there are special risk factors in the cochlear implant population, or whether different cochlear implant models pose different risks.

Adults and parents of children who are considering a cochlear implant or who have received cochlear implants should be advised of the risk of meningitis. They should also be informed of the availability of vaccines that have been shown to substantially reduce the incidence of meningitis in the general population resulting from the organisms that commonly cause bacterial meningitis (*Streptococcus pneumoniae*, *Haemophilus influenzae*, *Meningococcus*). National health agencies frequently provide updated information on the safety and utility of specific vaccines and offer recommendations reflecting local or regional conditions. Physicians or patients should refer to the applicable authorities for this information. These vaccines can be administered by pediatricians, primary care/family physicians, and infectious disease specialists.

Adults and parents of children who have received cochlear implants should be counseled on the symptoms of meningitis, the need to seek immediate medical care if any symptoms appear, and the need to advise the treating physicians of the presence of the cochlear implant and of the possibility of increased risk of meningitis associated with implant. They should also be counseled to obtain medical care at the first signs of otitis media.

- **Extreme direct pressure** on the implanted device, up, down, left or right may cause the implant to move and possibly dislodge the electrode array.
- **A direct impact to the implant site** may damage the implant and result in its failure to function. There have been instances of Advanced Bionics device failure as a result of a child hitting his/her head at the site of the implanted device. None of these reported incidents have resulted in a concussion or fracture of the skull. In all cases, the failed device was explanted and a new device reimplanted with no further complications.
- The long term effects of **chronic electrical stimulation** are unknown. Clinical experience with the system since 1991 has shown no adverse effects of chronic electrical stimulation on patient performance, electrical thresholds, or dynamic range.
- **Electrode displacement** can occur if the electrode is not inserted properly. Surgeons should be proficient in the use of the electrode insertion tool. **Failure to follow the recommended surgical procedure for placement and stabilization of the HiRes Ultra implant increases the risk of device migration or extrusion, and of damage resulting from impact trauma, including breakage of the electrode lead wires.**

Creating a recessed bed for the implant and securely stabilizing the device in place are critical elements of the surgical procedure.

- **Insertion of a cochlear implant electrode** will likely result in the loss of any residual hearing in the implanted ear.

In general, external components (e.g. sound processor and headpiece) should be removed or deactivated before undergoing the following medical procedures, where the implant may be exposed to electric fields or electrical current. If unexpected sounds or interruptions are experienced during or after medical procedures notify the physician.

- **Electrosurgery:** Electrosurgical instruments and RF Ablation instruments are capable of producing radio-frequency voltages of such magnitude that a direct coupling might occur between the cautery tip and the electrode. Induced currents could cause damage to the cochlear tissues or permanent damage to the implant. Monopolar electrosurgery and RF Ablation must only be used outside the head or neck region. For the use of bipolar electrosurgical equipment, the probe tips must not contact the implant and should be kept more than 1 mm (0.04 in) from the implant.
- **Extracorporeal Electrical Stimulation:** Electrodes or probe tips used in extracorporeal electrical stimulation instruments such as TENS, TEMS, etc. must not be placed over the implant.
- **Neurostimulation:** Do not use neurostimulation directly over the implant. High currents induced into the electrode lead can cause tissue damage or permanent damage to the implant.
- **Therapeutic Diathermy:** Therapeutic diathermy may induce currents in the electrode, which could cause injury to cochlear tissues. The HiRes Ultra 3D cochlear implant is safe to use during exposure to shortwave ($f=27.12$ MHz and $E\text{-field} = 1.16$ kV/m; $M\text{-field} = 8.5$ A/m) and longwave ($f = 1$ MHz and $E\text{-field} = 750$ V/m; $M\text{-field} = 2.0$ A/m) diathermy.
- **Ultrasound:** Diagnostic ultrasound can be used anywhere within the following parameters: Intensity (ISPTA): 1500 mW/cm², Frequency Range: 3.325 – 3.675 MHz, Duty Cycle: 20%. Therapeutic ultrasound energy must not be used over the implant since it may cause permanent damage to the implant by inadvertently concentrating the ultrasonic field and may cause injury in tissue surrounding the implant.
- **Electroconvulsive therapy:** Electroconvulsive therapy must never be used on a cochlear implant patient. Electroconvulsive therapy may cause tissue damage to the cochlea or permanent damage to the implant.
- **Ionizing Radiation Therapy:** The implant remained functional after being exposed to a total dosage of 250 Gray using a 15 MeV beam strength and 3 cm depth. Ionizing radiation therapy may damage the device. Long term effects of ionizing radiation on the implant may not be immediately detectable. The implant should not be placed directly in the ionizing radiation beam to minimize the risk of tissue necrosis. Diagnostic radiation imaging such as, CT, X-Ray, Mammography, etc. are safe to use.

MRI Warnings

Do not allow patients with a HiRes™ Ultra 3D cochlear implant to be in the area of an MRI scanner unless the following conditions have been met:

- The external sound processor and headpiece must be removed before entering a room containing an MRI scanner.
- Verify that the implant, or both implants if bilaterally implanted, are compatible for conducting an MRI before proceeding. Failure to do so can lead to device movement, device damage, magnet movement, patient discomfort, or trauma and pain to the patient.
- The recommended minimum duration of time post implant surgery prior to undergoing an MRI scan is 2 to 4 weeks in order to allow any inflammation to subside
- An MRI scan is not recommended if the patient has a fever.

Note: MRI procedures are contraindicated for CLARION (C1 and CII) cochlear implant recipients. For information regarding MRI use with HiRes™ 90K, HiRes™ 90K Advantage, and HiRes™ Ultra cochlear implants, please contact Advanced Bionics Technical Support.

Note: MRI safety was evaluated only for the HiRes™ Ultra 3D. Interactions between non-Advanced Bionics implants and the HiRes™ Ultra 3D during MRI are unknown.

Caution:

During the MRI procedure, you may experience pain, pressure or discomfort. If this occurs, please notify your physician.

Please consult with your physician prior to MRI to determine if the benefits of MRI are worthwhile over other imaging techniques.

MRI SAFETY INFORMATION

Testing has demonstrated that the HiRes™ Ultra 3D cochlear implant is MR Conditional. Unilateral and bilateral recipients with this device can be safely scanned in a horizontal closed bore quadrature coil MR system meeting the following conditions:

| MRI Field Strength | 1.5T | 3.0T |
|---------------------------------|-----------|----------|
| Maximum Spatial Field Gradient | 20 T/m | |
| RMS Gradient Field | 34.4 T/s | |
| Peak Slew Rate | 200 T/m/s | |
| Maximum whole body averaged SAR | 2.0 W/kg | 2.0 W/kg |
| Maximum head averaged SAR | 3.2 W/kg | 2.6 W/kg |

When tested under the scan conditions defined above the HiRes™ Ultra 3D cochlear implant produced a maximum temperature rise of <3°C after 15 minutes of continuous scanning.

Note: During the scan, patients might perceive auditory sensations. Adequate counseling of the patient is advised prior to performing the MRI. The likelihood and intensity of the auditory sensations can be reduced by selecting sequences with a lower Specific Absorption Rate (SAR) and slower gradient slew rates.

The largest artifact at 3.0T is >15 cm when imaged under spin echo and gradient echo sequences. Smaller artifacts are possible if using different scanning parameters or a non-magnetic plug.

Note: For cases that would clinically benefit from reduced device artifact (for example, some head or neck scans), the internal magnet is surgically removed and possibly replaced with a Temporary Non-Magnetic Plug before the recipient undergoes an MRI procedure. The HiRes™ Ultra 3D cochlear implant can withstand 5 replacement cycles.

For additional information regarding the use of an MRI scanner with a HiRes™ Ultra 3D cochlear implant, please contact Advanced Bionics Technical Support at technicalservices@advancedbionics.com or visit www.advancedbionics.com/mri.

PRECAUTIONS

- **Electrostatic Discharge (ESD):** It is known that static electricity can potentially damage sensitive electronic components such as the ones used in the cochlear implant system. Care should be taken to avoid situations in which high levels of static electricity are generated. More information is provided in the user manuals of the system. If static electricity is present, static electrical potential of the cochlear implant recipients can safely be reduced by the patients touching any person or object with their fingers prior to that person or object contacting the implant system.
- **Digital Cellular Phones:** Using or being in close vicinity to someone using some digital cellular phones may cause interference with the system. If such interference occurs, patients can turn off the sound processor or move a greater distance from the phone. Before purchasing a digital cellular phone, patients should evaluate whether it will interfere with their system. No such interference has been noted with cellular phones using analog technology.
- **Ingestion of Small Parts:** The external components of the implant system contain small parts that may be harmful if swallowed.
- **Airport/Security Metal Detectors:** Metal detectors, x-ray machines, and security scanners will not damage the implant or sound processor. However, individuals with a cochlear implant should be advised that passing through security metal detectors may activate the detector alarm. It is advised that patients carry their "Medical Device Identification Card" with them at all times. Cochlear implant users also might hear a distorted sound caused by the magnetic field around the security scanner door or hand-held scanning wand. Turning the sound-processor volume down before passing through security screening will ensure that those sounds, if they occur, are not too loud or uncomfortable.
- **Electromagnetic Interference:** RF workers may be exposed to higher interference. In the presence of high intensity EMI you may experience loss of sound. If this occurs, move from the area or temporarily discontinue use of the system by removing the headpiece.
- **Electromagnetic Sources:** It is advised to maintain a minimum distance of 30 cm (12 in) from electromagnetic emitters such as RFID, and metal detectors. It is advised to maintain a minimum distance of 50 cm (20 in) from an electronic article surveillance emitter.

- **Use of Another Person's Sound Processor:** Implant recipients should use only the sound processor that has been specifically programmed for them by their clinician. Use of a different sound processor may be ineffective in providing sound information and may cause physical discomfort.
- **Physical Activity:** When engaging in physical activities that include the possibility of trauma or impact, precautions should be taken, such as wearing a protective helmet, to reduce the risk of damage to the internal device.

EXPECTED LIFE

The expected life of the HiRes™ Ultra 3D cochlear implant is 10 years.

CLINICAL BENEFITS

The Intended Clinical Benefit of the HiRes™ Ultra 3D cochlear implant as part of the HiResolution™ Bionic Ear System is to:

- provide useful hearing to individuals with severe-to-profound hearing loss via electrical stimulation of the auditory nerve.
- provide a combined type of stimulation, electric plus acoustic, when sound processors are enabled with the Acoustic Earhook.

CLINICAL STUDIES

Performance Data

The HiRes™ Ultra 3D implant supports the HiResolution family of sound processing strategies including HiRes, HiRes with Fidelity 120 (HiRes 120), and ClearVoice.

HiRes and HiRes 120 Sound Processing

A clinical study was conducted in 50 adults implanted with a CII/HiRes 90K device who used a Harmony processor to document the benefits of HiRes 120 and HiRes sound processing. Performance with HiRes was assessed at the baseline visit and compared with HiRes 120 performance after three months of listening experience. Subsequently, subjects were refit and retested with HiRes. Results showed equivalent mean CNC word recognition scores for the two strategies. The mean HINT sentence perception scores in quiet and noise were significantly higher for HiRes 120 compared to baseline with HiRes. For HINT sentences in noise, the mean scores for HiRes 120 were significantly higher than scores after subjects were refit with HiRes.

Table 1. Mean Speech Scores for HiRes and HiRes 120

| Sound Processing Group | HiRes | HiRes 120 | HiRes |
|-------------------------------------|----------|-----------|----------|
| Test Interval | Baseline | 3 Months | 3 Months |
| CNC Words | 63 | 65 | 63 |
| HINT Sentences in Quiet | 88 | 93* | 91 |
| HINT Sentences in Noise (+8 dB SNR) | 64 | 70** | 65 |

* HiRes 120 score significantly different from baseline HiRes score ($p < .05$)

** HiRes 120 score significantly different from baseline and 3-month HiRes scores ($p < .05$)

Forty-three of 50 subjects (86%) preferred HiRes 120 over HiRes. Subjects rated strength of preference for the two strategies on a scale from 1 (weak preference) to 10 (strong preference). The mean strength of preference for the 43 subjects who preferred HiRes 120 was 7.9 (range: 1-10). The strength of preference was rated as 8 or higher by 26 of the 43 subjects, and 16 of the 43 subjects rated their preference as 10 (strong preference). For the 7 subjects who preferred HiRes, the mean strength of preference was 4.4 (range: 1-9).

ClearVoice

A clinical study was conducted in 46 adults who had at least six months experience with HiRes 120 sound processing and at least moderate speech perception abilities to investigate the benefits of ClearVoice. ClearVoice has three adaptive gain settings that allow individuals to select the setting that provides the best hearing—Low, Medium, and High. A two-week randomized, design was used to evaluate ClearVoice-Medium and ClearVoice-High. ClearVoice-Low was evaluated acutely during an initial test session. Speech benefit was compared for ClearVoice vs. HiRes 120 without ClearVoice (Control) using the AzBio sentence test.

Speech understanding in speech-spectrum noise was significantly better with ClearVoice-Medium and ClearVoice-High compared to the Control ($p < .0001$). ClearVoice-Medium significantly improved speech

understanding in multi-talker babble ($p < .02$). Speech understanding was no worse than the Control when listening in quiet for both ClearVoice-Medium and ClearVoice-High ($p < .0001$). Speech understanding with ClearVoice-Low was no worse than the Control in quiet, in speech-spectrum noise, and in multi-talker babble ($p < .001$).

Table 2. Mean AzBio Sentence Scores for HiRes 120 with and without ClearVoice

| Study Group | Control | ClearVoice Low | Control | ClearVoice Medium | Control | ClearVoice High |
|-----------------------|---------|----------------|---------|-------------------|---------|-----------------|
| Quiet | 87.3 | 87.8 | 88.6 | 88.3 | 86.8 | 87.7 |
| Speech-Spectrum Noise | 48.0 | 55.6 | 49.5 | 58.2 | 47.7 | 58.3 |
| Multi-Talker Babble | 42.8 | 47.2 | 44.9 | 48.1 | 44.9 | 46.2 |

Preference ratings indicated that 42 out of 45 subjects (93%) preferred ClearVoice to the Control for everyday listening (one subject did not complete the questionnaire). The mean strength of preference for the 42 subjects who preferred ClearVoice was 7.9 (1 = weak preference, 10 = strong preference).

Of the 42 subjects preferring ClearVoice, 22 indicated they would use it all of the time, 17 indicated they would use it most of the time, and 3 indicated they would use it some of the time. Of the 3 subjects preferring the Control, all indicated they would use ClearVoice some of the time.

NOTE: ClearVoice is only available in markets where ClearVoice has received regulatory approval. Contact Advanced Bionics for more information.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

A summary of safety and clinical performance (SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI: 08400944C160104Q2 and 08400944C160105Q4

<https://ec.europa.eu/tools/eudamed>

Users of the HiRes™ Ultra 3D cochlear implant in the European Union should report any serious incident to their local competent authority and to Advanced Bionics.

POSSIBLE ADVERSE EVENTS

The following risks associated with cochlear implantation and ear surgery also can occur.

- Implant patients incur the normal risks of surgery and general anesthesia.
- Major ear surgery may result in numbness, swelling or discomfort about the ear, disturbance of taste or balance, or neck pain. If these events occur, they are usually temporary and subside within a few weeks of surgery.
- Rarely, cochlear implantation may cause a leak of the inner ear fluid, which may result in meningitis.
- During the surgery, it is a rare possibility that the facial nerve could be injured resulting in a temporary or permanent weakening or full paralysis on the same side of the face as the implant.
- During the surgery, there is a rare possibility that cerebrospinal fluid leakage or perilymph fluid leakage could occur.
- As a result of the surgery, it is possible that dizziness, tinnitus, or vertigo may result. If these events occur, they are usually temporary and subside over time.
- The presence of a foreign body may cause irritation, inflammation, or skin breakdown and may require additional medical treatment or removal of the internal device.
- Skin infection in the area of the implant may require additional medical treatment or removal of the internal device.
- There is a possibility that the electrode or device may migrate requiring additional medical treatment or removal of the internal device to address any resulting injury.

PATIENT COUNSELING INFORMATION

Prospective cochlear implant candidates must be counseled appropriately on expected outcomes prior to surgery. Patients demonstrate a range of cochlear implant benefit.

Although it is not possible to predict post-implant performance preoperatively for individual patients, research and clinical experience have shown that age at implant, duration of severe-to-profound hearing loss, and preoperative speech perception skills have a significant effect on post-implant performance. Ear selection for implantation is left to the discretion of the patient, surgeon, and audiologist. There is no consensus in the field regarding implantation of the better versus poorer ear. If the poorer ear is implanted, patients should be counseled that postoperative performance ear may not equal that of the better non-implanted ear, especially if there also is long duration of deafness and negligible residual hearing preoperatively.

Communication mode (oral versus total communication) and the patient's auditory environment can affect outcomes in children. Implant-center professionals should counsel parents about the impact of communication mode and auditory environment on potential implant benefit in the pediatric population.

INFORMATION FOR USE AND REQUIRED TRAINING

A Surgeon's Manual and a video describing the surgical procedure and insertion of the electrode are provided to all ENT surgeons prior to implantation. ENT surgeons must be well versed in mastoid surgery and the facial recess approach to the round window. Advanced Bionics conducts periodic training courses on the recommended surgical procedure to implant HiRes™ Ultra. ENT surgeons who implant adults and/or pediatrics must receive training before implantation. ENT surgeons should work with an audiology professional who has been trained fully on the proper fitting and adjustment of the system. Failure to obtain the appropriate training might result in higher incidence of surgical and medical complications.

Device and Fitting Manuals are provided to all clinical centers with the Clinician's Programming System. Audiologists have experience in administering test procedures used to determine cochlear implant candidacy. They should be knowledgeable about state-of-the-art hearing systems and fitting procedures. In addition, at least one audiologist from a clinical center should be fully trained in the fitting of the Advanced Bionics cochlear implant in adults and/or children. Advanced Bionics conducts training courses for audiologists and strongly recommends that audiologists attend a training course. Failure to obtain the appropriate training might result in less-than-optimal patient performance.

Sound processor user guides are provided to all HiResolution™ Bionic Ear System recipients upon delivery of the system. Patient counseling materials are made available to all implant centers upon request. These materials provide detailed information about the system, indications for use, benefits, risks, and what is involved in patient selection, surgery and follow-up procedures.

SPECIAL HANDLING

Telemetry

The HiResolution™ Bionic Ear System incorporates bi-directional telemetry that verifies system function and continuously monitors the system during normal use.

Storage

The HiRes™ Ultra 3D implant should be stored at temperatures in the range of -20° to 33° Centigrade (4° to 91.4° Fahrenheit). Storage pressure range of 70 kPa to 150 kPa.

Handling

The HiRes™ Ultra implant package should be handled with care. An impact that damages the storage pack also could rupture the sterile packaging. The implant is not intended to tolerate a drop onto a hard surface. If the implant falls onto a hard surface, it must be returned to Advanced Bionics and the backup implant must be used. The returned device should be accompanied by a detailed description of the impact that caused the return.

Shelf Life

A Use Before date is noted on the packaging and is based on the date of the original sterilization.

Sterilization

The HiRes™ Ultra 3D cochlear implants are supplied in ethylene oxide sterile packaging with indicators of sterilization. Sterile packs should be inspected carefully to confirm that they have not been ruptured. Sterility cannot be guaranteed if the sterile package is damaged or opened. If the sterile pack of the HiRes™ Ultra 3D cochlear implant is damaged, the device must not be used and should be returned to Advanced Bionics. If the sterile packaging of the HiRes™ Ultra 3D cochlear implant has been opened, the cochlear implant cannot be resterilized by either the customer or Advanced Bionics.

Pressurized Environment Information

The HiRes Ultra 3D cochlear implants can withstand a pressure up to a depth of 42m under water (138 feet) or a gauge pressure of 4ATM (413 kPa).

Compatibility

The HiRes™ Ultra 3D cochlear implants are compatible with the following external system components and their corresponding accessories:

| | |
|-------------------------------------|--|
| HEADPIECES | Slim HP AquaMic™ Slim HP Slim HP Mic Slim HP Standard AquaMic™ 3D UHP 3D Plus UHP 3D |
| SOFTWARE | Target CI SoundWave™ Professional Suite Software (version 3.1 or later) |
| SOUND PROCESSORS¹ | Naïda™ CI M Sky CI™ M Naïda CI Q Neptune™ Sound Processor Harmony™ Processor |

¹ Requires CPI-3

NOTE: Compatible products are only available in markets that have received regulatory approval. Contact Advanced Bionics for more information.

| REF | PRODUCT NAME |
|------------|---|
| CI-1601-04 | HiRes™ Ultra 3D CI HiFocus™ MS Electrode |
| CI-1601-05 | HiRes™ Ultra 3D CI HiFocus™ SlimJ Electrode |
| CI-4509 | HiRes™ Ultra Reusable Surgical Tool Kit |
| CI-4331 | HiRes™ Ultra Recess Gauge |
| CI-4341 | HiRes™ Ultra Coil Gauge |
| CI-4421 | HiRes™ Ultra BTE Template |
| CI-1419 | HiRes™ Ultra 3D Replacement Magnet |
| CI-1420 | HiRes™ Ultra 3D Temporary Non-Magnetic Plug |
| CI-4349 | HiRes™ Ultra Magnet Pusher Tool |
| CI-4426 | HiRes™ Ultra Mock Up |
| CI-4508 | HiFocus™ Mid-Scala Electrode Instrument Kit |
| CI-4254 | HiFocus™ Mid-Scala Claw Tool |
| CI-4347 | HiFocus™ Mid-Scala Cochleostomy Gauge |
| CI-4507 | HiFocus™ MS Electrode Insertion Tool Kit |
| CI-4207 | HiFocus™ MS Electrode Insertion Tool |
| CI-1605 | HiFocus™ SlimJ Electrode Depth Gauge |
| CI-4350-02 | HiFocus™ Electrode Forceps Kit |

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