ADVANCED BIONICS

AB Remote App

L’application AB Remote

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CI-6058-001

instructions for use
Instructions d’utilisation
Instrucciones de uso
English, Français, Español

SYMBOLS TABLE:

Labeling symbols and their meanings:

	Date of Manufacture
	Manufacturer
	Model number
	Caution (refer to “Cautions & Warnings” section for further information)

LIMITATIONS, CONTRAINDICATIONS:

There are no known limitations or contraindications for the AB Remote app.

PRECAUTIONS, CAUTIONS, AND WARNINGS:

- *Do not allow children to play with or operate the app unattended.*
- *Do not attempt to pair your sound processor or hearing aids to any accessories while driving or operating heavy machinery.*
- *Do not stream audio to your sound processor or hearing aids while driving or operating heavy machinery.*
- *Use your sound processor for hands-free phone use only where permitted by law and only when you will not be distracted from the safe operation of your motor vehicle.*
- *Portable and mobile RF communications equipment, including radios and cellular phones, may affect sound quality of the Advanced Bionics Naída CI M and Sky CI M sound processors and accessories; however, there is no safety hazard associated with such equipment.*
- *Adjusting ambient balance while streaming audio (either with the multifunction button or with the mobile application) may interfere with the ability to be aware of surroundings.*

- *Do not attempt to service or modify the AB Remote App. Doing so may compromise system performance and will void the manufacturer's warranty. Products should be serviced only at Advanced Bionics and damaged products should be returned to Advanced Bionics.*

UNDESIRABLE SIDE EFFECTS:

No undesirable side effects.

PRODUCT DESCRIPTION AND PERFORMANCE CHARACTERISTICS:

The AB Remote is a smartphone app that is designed to provide user control over volume and program changes as well as provide access to help information.

The AB Remote app provides direct connectivity and control of Advanced Bionics Naída CI M and Sky CI M sound processors and Phonak Naída Link M and Sky Link M hearing instruments.

INTENDED USE AND USE ENVIRONMENT:

The AB Remote app is designed to be used with a Advanced Bionics Naída CI M and Sky CI M sound processors and Phonak Naída Link M and Sky Link M hearing instruments to allow the user to control volume related and program settings, as well as accessing additional helpful information. The intended use environments of the AB Remote app are healthcare and daily living environments.

EN

TARGET AUDIENCE/INTENDED USERS:

The intended users of the AB Remote App are recipients and/or care givers of Advanced Bionics Naída CI M and Sky CI M sound processors.

COMPATIBILITY:

The AB Remote app is compatible with the following products

- **Advanced Bionics Naída CI M and Sky CI M sound processors**
- **Phonak Naída Link M and Sky Link M hearing instruments**

The AB Remote app is compatible and has been tested with the following operating systems

- **Apple iOS 15 and later**
- **Android 8 and later**

The operating system compatibility information should be used as a reference guide only, because hardware and software levels and releases may vary, and numerous hardware and software combinations are possible. AB assumes no responsibility for the compatibility thus it is the user's responsibility to do product compatibility testing by himself before any purchasing decision is taken.

INSTALLATION INSTRUCTIONS:

Installation

Installing AB Remote is very simple. It installs just like other apps that are installed on your phone.

1. Turn on your smart phone.
2. Open your app store app or Google Play app.
3. Search for AB Remote or Advanced Bionics
4. Select the app and install it using your store account

USE INFORMATION:

Pairing

First use

The very first time you install the app, you will be guided to pair your devices. In 3 simple steps


1. Turn your devices off and then on to put them into pairing mode. View the instructions provided in the app if you are not sure how to do this
2. View your devices in the list of available devices
3. Tap “Connect” to pair your devices

**** NOTE: If you have a hearing device on the left and right side, they are paired and connected together in one step.***

Changing or Removing Paired Devices

You may have a need to change or remove the devices that are paired. For example, if you switch phones, or if you have replacement hearing devices.

1. Tap the icon in the top right corner for settings and help

2. Tap “My Devices”

3. Tap “Remove Devices”


4. Tap “Forget Devices” button shown in red

5. You may now pair new devices as described above.

Getting started with tips

The AB Remote app includes helpful animated tips.

Tips will be displayed the first time you run the app. At any time, you can access tips again by clicking this icon in the top left corner of the app:


Most of the common tasks you will need to do are explained in the tips section:

- **How to change the volume**
- **How to split the volume and control left and right side independently**
- **How to change the program**

- **How to control the balance between surrounding sounds and streaming sources like TV**
- **Additional Help**

Additional help, configuration, resources and information

Additional help, configuration, resources and information can be found by clicking the icon in the top right corner of the app.


From there, the following menu items are available.


- **“App settings”**
- **“My Devices”**
- **“Help” (link to your online account on My AB Online)** ***requires internet**
- **“Instructions For Use”** ***requires internet**

Device Check

The AB Remote app includes a screen to help you check the status of your devices.

To access the Device Check screen tap on the icon in the top right corner of AB Remote and tap “My Devices”

Speak at a normal level for 5 seconds at a distance of approximately 1-2 feet from the Sound Processor.

If the microphones are receiving sound, you will see the blue bars indicate audio and that your microphones are working correctly.

The Battery indicator will show the level of battery life remaining.

The Lock indicator will display the lock status of the processor.

This area will display blue bars when audio is detected by the microphones.

This area will display the status of the T-Mic.

The “Play Beep” button will play a beep that can be heard by the CI user

The “Blink Light” will cause the lights on the processor to illuminate

“Progressive Level”: If your healthcare professional has enabled the progressive levels feature, you can adjust the overall loudness in steps using the +/- buttons. The settings made here will continue even after powering off the cochlear implant.

Phonak Naída Link M and Sky Link M

If your device is Phonak Naída Link M and Sky Link M hearing Instrument, you will see indicators of the microphone status only.

Speak at a normal level for 5 seconds at a distance of approximately 1-2 feet from the Phonak Naída Link M and Sky Link M.

If the microphones are receiving sound, you will see the blue bars indicate audio

DISPOSAL:

The AB Remote app can be removed at any time using existing standard operating system features to remove any app. The AB Remote app does not collect or store patient data. No special provisions are required to remove patient data.

PERFORMANCE CHARACTERISTICS:

Not applicable.

CLINICAL BENEFITS:

N/A

PRODUCT OPTIONS:

Not applicable.

www.advancedbionics.com/clinicalandsafetyperformance/


Users of the AB Remote in the European Union should report any serious incident to their local competent authority.

GUIDELINES FOR IT SECURITY

Intended Use Environment

The AB Remote is intended for use in healthcare and home environments. The home environment is extended to include use outdoors and during travel (e.g., airplanes). The AB Remote application is intended to operate with Advanced Bionics Naída CI M and Sky CI M sound processors and Phonak Naída Link M and Sky Link M hearing instruments. The user is responsible for security of connections with other devices, including pairing the sound processor with other devices in a secure environment.

System and Network Diagram

System Interfaces & Connectivity

Component	Interface	Communication Protocol
Advanced Bionics Naída CI M sound processors **Advanced Bionics Sky CI M sound processors** **Phonak Naída Link M hearing instruments** **Phonak Sky Link M hearing instruments**	Bluetooth 4.2	Bluetooth 4.2

Cybersecurity Configuration

There is no specific AB Remote configuration needed to help ensure cybersecurity. Ensure pairing with Bluetooth devices is completed in a secure environment.

Integrated Cybersecurity Functions

Advanced Bionics Naída CI M and Sky CI M sound processors and AB Remote app do not include user authorization or authentication mechanisms. Each sound processor is uniquely configured to connect with a single Advanced Bionics cochlear implant as fitted by the Hearing Health Professional. Data integrity and quality of service is ensured through sound processor and implant error correction and detection mechanisms.

System, Configuration, and Data Backup/Restore

The AB Remote App does not provide System, Configuration or Data Backup/Restore functions. The user should visit his or her cochlear implant professional for any changes or restorations of sound processor settings.

Cybersecurity Incident Response

The AB Remote App does not detect cybersecurity attacks or intrusions. In case of abnormal behavior of the AB Remote app or sound processor, please contact Advanced Bionics.

Patches and Updates

Advanced Bionics may release software updates from time to time to address defects or enhance features. The software updates will be notified through the existing app store mechanisms.

Available Training

Training is not available for the AB Remote App.

End of Cybersecurity Support

Advanced Bionics Apps and Sound Processors are supported until they no longer have regulatory approval in the country of use or until the product is obsolesced by Advanced Bionics.

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©2025 Advanced Bionics LLC. All Rights Reserved /
No hay capacitación disponible para la aplicación AB Remote.
Advanced Bionics LLC