



Off-the-Ear Power Options

by Advanced Bionics®



Flexible. Comfortable. Secure. Convenient.



Advanced Bionics introduces off-the-ear power options for the Harmony® and Auria® BTE sound processors. **Small and lightweight**, the new off-the-ear choices give every recipient **flexibility and comfort, security and convenience**.



A lighter weight alternative

From active adults to tiny toddlers, less on-the-ear translates to more comfort for every recipient.

Peace of mind

Designed for active lifestyles, you can feel confident that the battery will stay securely in place.

Grows with your child

The new off-the-ear power options allow the Harmony system to grow with your child from infancy to adulthood. Start with a secure and lightweight off-the-ear option and transition to a full behind-the-ear system as your child grows.

Options for every adult


Whether you are an active adult wanting less on-the-ear or a bilateral recipient looking for a convenient power option, the new off-the-ear power options provide every adult the flexibility to match their lifestyle.

Bilateral Power*

Two processors, one battery. Only AB gives bilateral recipients of all ages the option to power both processors with a single PowerCel™ adapter.


* The PowerCel Adapter may not have sufficient power to support use by all bilateral implant recipients.

Off-the-Ear Power Option



- Harmony BTE sound processor with BTE Power Adapter
- Cable Options: 5", 11", 22", 32"
- PowerCel Adapter with PowerCel attached.**

Bilateral Power Option*



- Two Harmony BTE sound processors with Bilateral Power Adapter and PowerCel attached**
- Bilateral Power Adapter
- PowerCel Adapter with PowerCel attached.**

**Compatible with Harmony and Auria sound processors.

HEADQUARTERS

Advanced Bionics, LLC
12740 San Fernando Road
Sylmar, CA 91342 USA
877. 829.0026 in US and Canada
800. 678.3575 TTY
661. 362.1400
661. 362.1500 Fax
info@AdvancedBionics.com

EUROPE

Advanced Bionics SARL
76 rue de Battenheim
68170 Rixheim, France
+33.0.3.89.65.98.00
+33.0.3.89.65.50.05 Fax
europe@AdvancedBionics.com

ASIA PACIFIC

Advanced Bionics Asia-Pacific Limited
Suite 4203, 42/F, Tower One
Lippo Centre, 89 Queensway
Hong Kong
852.2526.7668
852.2526.7628 Fax
AP@AdvancedBionics.com

LATIN AMERICA

Advanced Bionics
28515 Westinghouse Place
Valencia, CA 91355, USA
661. 362.1840
661. 362.4604 Fax
LA@AdvancedBionics.com

For Insurance Reimbursement Assistance:

AB Insurance Reimbursement Specialists are available to assist you in determining your health plans coverage. If your plan requires pre-authorization, AB Insurance Reimbursement Specialists will coordinate with your provider to obtain the information required for approval. AB cannot guarantee coverage or payment by your health plan.

Email: insurance@AdvancedBionics.com
Phone: 877.779.0229
TTY: 1.800.678.3575
Fax: 877.833.6318



To place your order:

Online: www.BionicEar.com/webstore
Email: customerservice@AdvancedBionics.com
Webchat: www.BionicEar.com/support
Phone: 877.829.0026
TTY 800.678.3575